

# Product Datasheet

## SPOP Antibody (S01-5B7)

### NBP3-19781-100ul

Unit Size: 100 ul

Store at 4C short term. Aliquot and store at -20C long term. Avoid freeze-thaw cycles.

[www.novusbio.com](http://www.novusbio.com)



[technical@novusbio.com](mailto:technical@novusbio.com)

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:  
[www.novusbio.com/NBP3-19781](http://www.novusbio.com/NBP3-19781)

Updated 4/15/2024 v.20.1

Earn rewards for product  
reviews and publications.

Submit a publication at [www.novusbio.com/publications](http://www.novusbio.com/publications)

Submit a review at [www.novusbio.com/reviews/destination/NBP3-19781](http://www.novusbio.com/reviews/destination/NBP3-19781)



**NBP3-19781-100ul**

SPOP Antibody (S01-5B7)


| Product Information     |                                                                                        |
|-------------------------|----------------------------------------------------------------------------------------|
| Unit Size               | 100 ul                                                                                 |
| Concentration           | 0.3 mg/ml                                                                              |
| Storage                 | Store at 4C short term. Aliquot and store at -20C long term. Avoid freeze-thaw cycles. |
| Clonality               | Monoclonal                                                                             |
| Clone                   | S01-5B7                                                                                |
| Preservative            | 0.01% Sodium Azide                                                                     |
| Isotype                 | IgG                                                                                    |
| Purity                  | Affinity purified                                                                      |
| Buffer                  | 50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.05% BSA                         |
| Target Molecular Weight | 42 kDa                                                                                 |

| Product Description |                                                      |
|---------------------|------------------------------------------------------|
| Host                | Rabbit                                               |
| Gene Symbol         | SPOP                                                 |
| Species             | Human                                                |
| Immunogen           | A synthetic peptide of human SPOP (Uniprot # O43791) |

| Product Application Details |              |
|-----------------------------|--------------|
| Applications                | Western Blot |
| Recommended Dilutions       | Western Blot |

**Images**

Western Blot: SPOP Antibody (S01-5B7) [NBP3-19781] - Western blot detection of SPOP in MCF-7 using NBP3-19781. Predicted band size: 42kDa. Observed band size: 42kDa.





### **Novus Biologicals USA**

10730 E. Briarwood Avenue  
Centennial, CO 80112  
USA  
Phone: 303.730.1950  
Toll Free: 1.888.506.6887  
Fax: 303.730.1966  
nb-customerservice@bio-techne.com

### **Bio-Techne Canada**

21 Canmotor Ave  
Toronto, ON M8Z 4E6  
Canada  
Phone: 905.827.6400  
Toll Free: 855.668.8722  
Fax: 905.827.6402  
canada.inquires@bio-techne.com

### **Bio-Techne Ltd**

19 Barton Lane  
Abingdon Science Park  
Abingdon, OX14 3NB, United Kingdom  
Phone: (44) (0) 1235 529449  
Free Phone: 0800 37 34 15  
Fax: (44) (0) 1235 533420  
info.EMEA@bio-techne.com

### **General Contact Information**

[www.novusbio.com](http://www.novusbio.com)  
Technical Support: nb-technical@bio-techne.com  
Orders: nb-customerservice@bio-techne.com  
General: novus@novusbio.com

### **Products Related to NBP3-19781-100ul**

---

|                   |                                                     |
|-------------------|-----------------------------------------------------|
| HAF008            | Goat anti-Rabbit IgG Secondary Antibody [HRP]       |
| NB7160            | Goat anti-Rabbit IgG (H+L) Secondary Antibody [HRP] |
| NBP2-24891        | Rabbit IgG Isotype Control                          |
| NBP1-51077-0.05mg | Recombinant Human SPOP His Protein                  |

---

### **Limitations**

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit [www.novusbio.com/guarantee](http://www.novusbio.com/guarantee)

Earn gift cards/discounts by submitting a review: [www.novusbio.com/reviews/submit/NBP3-19781](http://www.novusbio.com/reviews/submit/NBP3-19781)

Earn gift cards/discounts by submitting a publication using this product:  
[www.novusbio.com/publications](http://www.novusbio.com/publications)

