

Product Datasheet

WDR1 Antibody (S01-6D6) NBP3-19802-100ul

Unit Size: 100 ul

Store at 4C short term. Aliquot and store at -20C long term. Avoid freeze-thaw cycles.

www.novusbio.com



technical@novusbio.com

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/NBP3-19802

Updated 9/9/2025 v.20.1

Earn rewards for product
reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/NBP3-19802



NBP3-19802-100ul


WDR1 Antibody (S01-6D6)

Product Information	
Unit Size	100 ul
Concentration	0.3 mg/ml
Storage	Store at 4C short term. Aliquot and store at -20C long term. Avoid freeze-thaw cycles.
Clonality	Monoclonal
Clone	S01-6D6
Preservative	0.01% Sodium Azide
Isotype	IgG
Purity	Affinity purified
Buffer	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.05% BSA
Target Molecular Weight	66 kDa
Product Description	
Description	Novus Biologicals Rabbit WDR1 Antibody (S01-6D6) (NBP3-19802) is a recombinant monoclonal antibody validated for use in IHC, WB and ICC/IF. All Novus Biologicals antibodies are covered by our 100% guarantee.
Host	Rabbit
Gene ID	9948
Gene Symbol	WDR1
Species	Human, Mouse, Rat
Immunogen	A synthetic peptide of human WDR1 (Uniprot # O75083)
Product Application Details	
Applications	Western Blot, Immunohistochemistry-Paraffin, Immunocytochemistry/Immunofluorescence, Immunohistochemistry
Recommended Dilutions	Western Blot 1/1000-1/5000, Immunohistochemistry 1/20-1/50, Immunocytochemistry/ Immunofluorescence 1/20-1/50, Immunohistochemistry-Paraffin




Images


Western Blot: WDR1 Antibody (S01-6D6) [NBP3-19802] - Western blot detection of WDR1 in Jurkat cell lysates using NBP3-19802 (1:500 diluted). Predicted band size: 66kDa. Observed band size: 66kDa.



Immunocytochemistry/Immunofluorescence: WDR1 Antibody (S01-6D6) [NBP3-19802] - Immunocytochemistry of WDR1 (green) in hela using NBP3-19802 at dilution 1/5, and DAPI (blue)



Immunohistochemistry-Paraffin: WDR1 Antibody (S01-6D6) [NBP3-19802] - Immunohistochemistry of WDR1 in paraffin-embedded Human Cholangiocarcinoma using NBP3-19802 at dilution 1/50



**Novus Biologicals USA**

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: nb-technical@bio-techne.com
Orders: nb-customerservice@bio-techne.com
General: novus@novusbio.com

Products Related to NBP3-19802-100ul

HAF008	Goat anti-Rabbit IgG Secondary Antibody [HRP]
NB7160	Goat anti-Rabbit IgG (H+L) Secondary Antibody [HRP]
NBP2-24891	Rabbit IgG Isotype Control
H00009948-P01-10ug	Recombinant Human WDR1 GST (N-Term) Protein

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/NBP3-19802

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

