

Product Datasheet

Histone H4 Antibody (S05-8H9) NBP3-20040-100ul

Unit Size: 100 ul

Store at 4C short term. Aliquot and store at -20C long term. Avoid freeze-thaw cycles.

www.novusbio.com



technical@novusbio.com

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/NBP3-20040

Updated 2/20/2025 v.20.1

Earn rewards for product
reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/NBP3-20040



NBP3-20040-100ul

Histone H4 Antibody (S05-8H9)

Product Information	
Unit Size	100 ul
Concentration	0.3 mg/ml
Storage	Store at 4C short term. Aliquot and store at -20C long term. Avoid freeze-thaw cycles.
Clonality	Monoclonal
Clone	S05-8H9
Preservative	0.02% Sodium Azide
Isotype	IgG
Purity	Affinity purified
Buffer	10mM PBS (pH 7.4), 150mM sodium chloride, 0.05% BSA, 50% glycerol
Target Molecular Weight	11 kDa


Product Description	
Host	Rabbit
Gene ID	121504
Gene Symbol	H4C16
Species	Human, Mouse, Rat
Immunogen	A synthetic peptide of human Histone H4 (Uniprot # P62805)

Product Application Details	
Applications	Western Blot, Chromatin Immunoprecipitation, Immunocytochemistry/ Immunofluorescence, Immunohistochemistry, Immunohistochemistry-Paraffin, Immunoprecipitation
Recommended Dilutions	Western Blot 1:1000-1:2000, Chromatin Immunoprecipitation 1:20-1:50, Immunohistochemistry 1:100-1:200, Immunocytochemistry/ Immunofluorescence 1:50-1:200, Immunoprecipitation 1:20-1:50, Immunohistochemistry-Paraffin




Images

Western Blot: Histone H4 Antibody (S05-8H9) [NBP3-20040] - Western blot detection of Histone H4 in Hela cell lysates using NBP3-20040 (1:1000 diluted). Predicted band size: 11kDa. Observed band size: 11kDa.



Immunohistochemistry-Paraffin: Histone H4 Antibody (S05-8H9) [NBP3-20040] - Immunohistochemistry of Histone H4 in paraffin-embedded Human colon cancer tissue using NBP3-20040 at dilution 1/100





Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: nb-technical@bio-techne.com
Orders: nb-customerservice@bio-techne.com
General: novus@novusbio.com

Products Related to NBP3-20040-100ul

HAF008	Goat anti-Rabbit IgG Secondary Antibody [HRP]
NB7160	Goat anti-Rabbit IgG (H+L) Secondary Antibody [HRP]
NBP2-24891	Rabbit IgG Isotype Control
NB21-2011PEP	Histone H4 [Monomethyl Arg3] Antibody Blocking Peptide

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/NBP3-20040

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

