

Product Datasheet

dsRNA Antibody (1D3) - Chimeric - Azide and BSA Free NBP3-20111-0.2mg

Unit Size: 0.2 mg

Store at 4C for up to 3 months. For longer storage, aliquot and store at -20C.

www.novusbio.com



technical@novusbio.com

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/NBP3-20111

Updated 9/9/2025 v.20.1

Earn rewards for product
reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/NBP3-20111



NBP3-20111-0.2mg

dsRNA Antibody (1D3) - Chimeric - Azide and BSA Free

Product Information	
Unit Size	0.2 mg
Concentration	1 mg/ml
Storage	Store at 4C for up to 3 months. For longer storage, aliquot and store at -20C.
Clonality	Monoclonal
Clone	1D3
Preservative	0.02% Proclin 300
Isotype	IgG Kappa
Purity	Protein A purified
Buffer	PBS


Product Description	
Description	Novus Biologicals Rabbit dsRNA Antibody (1D3) - Chimeric - Azide and BSA Free (NBP3-20111) is a recombinant monoclonal antibody validated for use in IHC, WB and ELISA. All Novus Biologicals antibodies are covered by our 100% guarantee.
Host	Rabbit
Species	Multi-Species
Immunogen	RDV-RNA-methylated bovine serum albumin complex.

Product Application Details	
Applications	Western Blot, ELISA, Immunohistochemistry
Recommended Dilutions	Western Blot, ELISA, Immunohistochemistry
Application Notes	This chimeric rabbit antibody was made using the variable domain sequences of the original Mouse format, for improved compatibility with existing reagents, assays and techniques. The antibody binds specifically to dsRNA, which is part of the genetic material of some viruses, such as rotaviruses and the bluetongue virus.




Images

Immunocytochemistry/ Immunofluorescence: Rabbit dsRNA mAb (1D3) - Chimeric - Azide and BSA Free [NBP3-20111] - Immunofluorescence analysis of Zika virus-infected A549 cells stained with the chimeric rabbit IgG version of 1D3 followed by a secondary antibody. The nuclear stain is DAPI (blue). Panels show from left-right: (green) with DAPI and merged channels (green), DAPI (blue), Zika's E protein (red)). Images are courtesy of Cecilia Alejandra Vazquez (IQUIBICEN (UBA-CONICET)).



Immunocytochemistry/ Immunofluorescence: Rabbit dsRNA mAb (1D3) - Chimeric - Azide and BSA Free [NBP3-20111] - Immunofluorescence analysis of paraformaldehyde-fixed, triton-permeabilized cells infected with Dengue 2 virus stained with the chimeric rabbit IgG version of 1D3 (Ab00460-23.0) at 1/300 followed by Alexa 568 secondary antibody (red). The image is courtesy of Cecilia Alejandra Vazquez (IQUIBICEN (UBA-CONICET)).





Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: nb-technical@bio-techne.com
Orders: nb-customerservice@bio-techne.com
General: novus@novusbio.com

Products Related to NBP3-20111-0.2mg

HAF008	Goat anti-Rabbit IgG Secondary Antibody [HRP]
NB7160	Goat anti-Rabbit IgG (H+L) Secondary Antibody [HRP]
NBP3-11403	dsRNA Antibody Pack

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/NBP3-20111

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

