

# Product Datasheet

## PRMT5 Antibody (RAB-C136) [PE] - Chimeric NBP3-20131PE

Unit Size: 0.1 ml

Store at 4C in the dark.

[www.novusbio.com](http://www.novusbio.com)



[technical@novusbio.com](mailto:technical@novusbio.com)

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:  
[www.novusbio.com/NBP3-20131PE](http://www.novusbio.com/NBP3-20131PE)

Updated 11/11/2025 v.20.1

Earn rewards for product  
reviews and publications.

Submit a publication at [www.novusbio.com/publications](http://www.novusbio.com/publications)

Submit a review at [www.novusbio.com/reviews/destination/NBP3-20131PE](http://www.novusbio.com/reviews/destination/NBP3-20131PE)



**NBP3-20131PE**

PRMT5 Antibody (RAB-C136) [PE] - Chimeric

| Product Information |                                                                                             |
|---------------------|---------------------------------------------------------------------------------------------|
| Unit Size           | 0.1 ml                                                                                      |
| Concentration       | Please see the vial label for concentration. If unlisted please contact technical services. |
| Storage             | Store at 4C in the dark.                                                                    |
| Clonality           | Monoclonal                                                                                  |
| Clone               | RAB-C136                                                                                    |
| Preservative        | 0.05% Sodium Azide                                                                          |
| Isotype             | IgG Kappa                                                                                   |
| Conjugate           | PE                                                                                          |
| Purity              | Protein A purified                                                                          |
| Buffer              | PBS                                                                                         |

| Product Description |                                                                                                                                   |
|---------------------|-----------------------------------------------------------------------------------------------------------------------------------|
| Host                | Rabbit                                                                                                                            |
| Gene ID             | 10419                                                                                                                             |
| Gene Symbol         | PRMT5                                                                                                                             |
| Species             | Human                                                                                                                             |
| Immunogen           | This antibody was obtained by recombinant antibody (rAb) phage display recognizing PRMT5 protein under non-denaturing conditions. |

| Product Application Details |                                                                                                                                              |
|-----------------------------|----------------------------------------------------------------------------------------------------------------------------------------------|
| Applications                | ELISA, Flow Cytometry, Immunocytochemistry/ Immunofluorescence, Immunoprecipitation, Chromatin Immunoprecipitation (ChIP), Mass Spectrometry |
| Recommended Dilutions       | Flow Cytometry, ELISA, Immunocytochemistry/ Immunofluorescence, Immunoprecipitation, Mass Spectrometry, Chromatin Immunoprecipitation (ChIP) |
| Application Notes           | Optimal dilution of this antibody should be experimentally determined.                                                                       |





### **Novus Biologicals USA**

10730 E. Briarwood Avenue  
Centennial, CO 80112  
USA  
Phone: 303.730.1950  
Toll Free: 1.888.506.6887  
Fax: 303.730.1966  
nb-customerservice@bio-techne.com

### **Bio-Techne Canada**

21 Canmotor Ave  
Toronto, ON M8Z 4E6  
Canada  
Phone: 905.827.6400  
Toll Free: 855.668.8722  
Fax: 905.827.6402  
canada.inquires@bio-techne.com

### **Bio-Techne Ltd**

19 Barton Lane  
Abingdon Science Park  
Abingdon, OX14 3NB, United Kingdom  
Phone: (44) (0) 1235 529449  
Free Phone: 0800 37 34 15  
Fax: (44) (0) 1235 533420  
info.EMEA@bio-techne.com

### **General Contact Information**

[www.novusbio.com](http://www.novusbio.com)  
Technical Support: [nb-technical@bio-techne.com](mailto:nb-technical@bio-techne.com)  
Orders: [nb-customerservice@bio-techne.com](mailto:nb-customerservice@bio-techne.com)  
General: [novus@novusbio.com](mailto:novus@novusbio.com)

### **Products Related to NBP3-20131PE**

---

|                    |                                              |
|--------------------|----------------------------------------------|
| H00010419-P01-10ug | Recombinant Human PRMT5 GST (N-Term) Protein |
| NBL1-11567         | H4/j Overexpression Lysate                   |
| H00010419-Q01-10ug | Recombinant Human PRMT5 GST (N-Term) Protein |
| NB600-235          | BLIMP1/PRDM1 Antibody (3H2-E8) - BSA Free    |

---

### **Limitations**

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit [www.novusbio.com/guarantee](http://www.novusbio.com/guarantee)

Earn gift cards/discounts by submitting a review: [www.novusbio.com/reviews/submit/NBP3-20131PE](http://www.novusbio.com/reviews/submit/NBP3-20131PE)

Earn gift cards/discounts by submitting a publication using this product:  
[www.novusbio.com/publications](http://www.novusbio.com/publications)

