

Product Datasheet

L3MBTL1 Antibody (RAB-C213) [FITC] - Chimeric NBP3-20137F

Unit Size: 0.1 ml

Store at 4C in the dark.

www.novusbio.com



technical@novusbio.com

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/NBP3-20137F

Updated 11/12/2023 v.20.1

Earn rewards for product
reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/NBP3-20137F



NBP3-20137F

L3MBTL1 Antibody (RAB-C213) [FITC] - Chimeric

| Product Information | |
|---------------------|---------------------------------------------------------------------------------------------|
| Unit Size | 0.1 ml |
| Concentration | Please see the vial label for concentration. If unlisted please contact technical services. |
| Storage | Store at 4C in the dark. |
| Clonality | Monoclonal |
| Clone | RAB-C213 |
| Preservative | 0.05% Sodium Azide |
| Isotype | IgG Kappa |
| Conjugate | FITC |
| Purity | Protein A purified |
| Buffer | PBS |

| Product Description | |
|---------------------|-------------------------------------------------------------------------------------------------------------------------------------|
| Host | Rabbit |
| Gene ID | 26013 |
| Gene Symbol | L3MBTL1 |
| Species | Human |
| Immunogen | This antibody was obtained by recombinant antibody (rAb) phage display recognizing L3MBTL1 protein under non-denaturing conditions. |

| Product Application Details | |
|-----------------------------|---------------------------------------------------------------------------------------------------------------------------------------------|
| Applications | ELISA, Flow Cytometry, Immunocytochemistry/Immunofluorescence, Immunoprecipitation, Chromatin Immunoprecipitation (ChIP), Mass Spectrometry |
| Recommended Dilutions | Flow Cytometry, ELISA, Immunocytochemistry/Immunofluorescence, Immunoprecipitation, Mass Spectrometry, Chromatin Immunoprecipitation (ChIP) |
| Application Notes | Optimal dilution of this antibody should be experimentally determined. |





Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: nb-technical@bio-techne.com
Orders: nb-customerservice@bio-techne.com
General: novus@novusbio.com

Products Related to NBP3-20137F

| | |
|--------------------|------------------------------------------------|
| H00026013-P01-10ug | Recombinant Human L3MBTL1 GST (N-Term) Protein |
| H00008362-P01-2ug | Recombinant Human H4/d GST (N-Term) Protein |
| NBP2-07169 | L3MBTL1 Overexpression Lysate |
| NBP1-96140 | BMI-1 Antibody (LLBmi1-1) - BSA Free |

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/NBP3-20137F

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

