

# Product Datasheet

## Fas/TNFRSF6/CD95 Antibody (rB-R18) [DyLight 550] NBP3-20704R

Unit Size: 0.1 ml

Store at 4C in the dark.

[www.novusbio.com](http://www.novusbio.com)



[technical@novusbio.com](mailto:technical@novusbio.com)

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:  
[www.novusbio.com/NBP3-20704R](http://www.novusbio.com/NBP3-20704R)

Updated 10/26/2023 v.20.1

Earn rewards for product  
reviews and publications.

Submit a publication at [www.novusbio.com/publications](http://www.novusbio.com/publications)

Submit a review at [www.novusbio.com/reviews/destination/NBP3-20704R](http://www.novusbio.com/reviews/destination/NBP3-20704R)



**NBP3-20704R**

Fas/TNFRSF6/CD95 Antibody (rB-R18) [DyLight 550]

| Product Information |   |
|---------------------|---|
| Unit Size           | 0.1 ml  |
| Concentration       | Please see the vial label for concentration. If unlisted please contact technical services. |
| Storage             | Store at 4C in the dark.  |
| Clonality           | Monoclonal  |
| Clone               | rB-R18  |
| Preservative        | 0.05% Sodium Azide  |
| Isotype             | IgG2a Kappa   |
| Conjugate           | DyLight 550   |
| Purity              | Protein A or G purified   |
| Buffer              | 50mM Sodium Borate  |

| Product Description     |  |
|-------------------------|--|
| Host                    | Mouse  |
| Gene ID                 | 355  |
| Gene Symbol             | FAS  |
| Species                 | Human  |
| Specificity/Sensitivity | monoclonal antibody rB-R18 specifically recognizes Fas/TNFRSF6/CD95, also known as Fas, a transmembrane glycoprotein with a MW of 40-45kDa, containing 8kDa of N-glycoside-linked polysaccharide |
| Immunogen               | Recombinant human Fas/TNFRSF6/CD95 protein   |
| Notes                   | DyLight (R) is a trademark of Thermo Fisher Scientific Inc. and its subsidiaries.  |

| Product Application Details |   |
|-----------------------------|---|
| Applications                | Flow Cytometry, Immunocytochemistry/Immunofluorescence, Immunohistochemistry-Frozen |
| Recommended Dilutions       | Flow Cytometry, Immunocytochemistry/Immunofluorescence, Immunohistochemistry-Frozen |
| Application Notes           | Optimal dilution of this antibody should be experimentally determined.              |





### **Novus Biologicals USA**

10730 E. Briarwood Avenue  
Centennial, CO 80112  
USA

Phone: 303.730.1950  
Toll Free: 1.888.506.6887  
Fax: 303.730.1966  
novus@novusbio.com

### **Bio-Techne Canada**

21 Canmotor Ave  
Toronto, ON M8Z 4E6  
Canada

Phone: 905.827.6400  
Toll Free: 855.668.8722  
Fax: 905.827.6402  
canada.inquires@bio-techne.com

### **Bio-Techne Ltd**

19 Barton Lane  
Abingdon Science Park  
Abingdon, OX14 3NB, United Kingdom

Phone: (44) (0) 1235 529449  
Free Phone: 0800 37 34 15  
Fax: (44) (0) 1235 533420  
info.EMEA@bio-techne.com

### **General Contact Information**

www.novusbio.com  
Technical Support: technical@novusbio.com  
Orders: orders@novusbio.com  
General: novus@novusbio.com

### **Products Related to NBP3-20704R**

---

|                |  |
|----------------|--|
| NBP1-96981R    | Mouse IgG2a Kappa Isotype Control (M2AK) [DyLight 550] |
| NBP2-61594-5ug | Recombinant Human Fas/TNFRSF6/CD95 Protein             |
| 210-TA-005     | TNF-alpha [Unconjugated]                               |
| 7398-FS-050    | Fas/TNFRSF6/CD95                                       |

---

### **Limitations**

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit [www.novusbio.com/guarantee](http://www.novusbio.com/guarantee)

Earn gift cards/discounts by submitting a review: [www.novusbio.com/reviews/submit/NBP3-20704R](http://www.novusbio.com/reviews/submit/NBP3-20704R)

Earn gift cards/discounts by submitting a publication using this product:  
[www.novusbio.com/publications](http://www.novusbio.com/publications)

