

# Product Datasheet

## Cytokeratin 18 Antibody (rKRT18/7299) [Biotin] NBP3-20911B

Unit Size: 0.1 ml

Store at 4C in the dark.

[www.novusbio.com](http://www.novusbio.com)



[technical@novusbio.com](mailto:technical@novusbio.com)

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:  
[www.novusbio.com/NBP3-20911B](http://www.novusbio.com/NBP3-20911B)

Updated 10/26/2023 v.20.1

Earn rewards for product  
reviews and publications.

Submit a publication at [www.novusbio.com/publications](http://www.novusbio.com/publications)

Submit a review at [www.novusbio.com/reviews/destination/NBP3-20911B](http://www.novusbio.com/reviews/destination/NBP3-20911B)



**NBP3-20911B**

Cytokeratin 18 Antibody (rKRT18/7299) [Biotin]

| Product Information         |   |
|-----------------------------|---|
| Unit Size                   | 0.1 ml  |
| Concentration               | Please see the vial label for concentration. If unlisted please contact technical services.   |
| Storage                     | Store at 4C in the dark.  |
| Clonality                   | Monoclonal  |
| Clone                       | rKRT18/7299   |
| Preservative                | 0.05% Sodium Azide  |
| Isotype                     | IgG1 Kappa  |
| Conjugate                   | Biotin  |
| Purity                      | Protein A or G purified   |
| Buffer                      | PBS   |
| Product Description         |   |
| Host                        | Mouse   |
| Gene ID                     | 3875  |
| Gene Symbol                 | KRT18   |
| Species                     | Human   |
| Specificity/Sensitivity     | This monoclonal antibody reacts with a wide variety of simple epithelia It does not react with stratified squamous epithelia It reacts with epithelial tumors of the gastrointestinal tract, lung, breast, pancreas, ovary, and thyroid |
| Immunogen                   | Recombinant full-length human Cytokeratin 18 protein  |
| Product Application Details |   |
| Applications                | Western Blot, Flow Cytometry, Immunocytochemistry/Immunofluorescence, Immunohistochemistry-Paraffin   |
| Recommended Dilutions       | Western Blot, Flow Cytometry, Immunocytochemistry/Immunofluorescence, Immunohistochemistry-Paraffin   |
| Application Notes           | Optimal dilution of this antibody should be experimentally determined.  |





### **Novus Biologicals USA**

10730 E. Briarwood Avenue  
Centennial, CO 80112  
USA  
Phone: 303.730.1950  
Toll Free: 1.888.506.6887  
Fax: 303.730.1966  
novus@novusbio.com

### **Bio-Techne Canada**

21 Canmotor Ave  
Toronto, ON M8Z 4E6  
Canada  
Phone: 905.827.6400  
Toll Free: 855.668.8722  
Fax: 905.827.6402  
canada.inquires@bio-techne.com

### **Bio-Techne Ltd**

19 Barton Lane  
Abingdon Science Park  
Abingdon, OX14 3NB, United Kingdom  
Phone: (44) (0) 1235 529449  
Free Phone: 0800 37 34 15  
Fax: (44) (0) 1235 533420  
info.EMEA@bio-techne.com

### **General Contact Information**

www.novusbio.com  
Technical Support: technical@novusbio.com  
Orders: orders@novusbio.com  
General: novus@novusbio.com

### **Products Related to NBP3-20911B**

---

|                    |   |
|--------------------|---|
| IC002B             | Mouse IgG1 Isotype Control (11711) [Biotin]           |
| H00003875-P01-10ug | Recombinant Human Cytokeratin 18 GST (N-Term) Protein |
| D6050              | IL-6 [HRP]  |
| H00003875-P02-10ug | Recombinant Human Cytokeratin 18 GST (N-Term) Protein |

---

### **Limitations**

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit [www.novusbio.com/guarantee](http://www.novusbio.com/guarantee)

Earn gift cards/discounts by submitting a review: [www.novusbio.com/reviews/submit/NBP3-20911B](http://www.novusbio.com/reviews/submit/NBP3-20911B)

Earn gift cards/discounts by submitting a publication using this product:  
[www.novusbio.com/publications](http://www.novusbio.com/publications)

