

Product Datasheet

Synaptophysin Antibody (SYP/7976R) [DyLight 650] NBP3-20998C

Unit Size: 0.1 ml

Store at 4C in the dark.

www.novusbio.com



technical@novusbio.com

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/NBP3-20998C

Updated 10/26/2023 v.20.1

Earn rewards for product
reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/NBP3-20998C



NBP3-20998C

Synaptophysin Antibody (SYP/7976R) [DyLight 650]

Product Information	
Unit Size	0.1 ml
Concentration	Please see the vial label for concentration. If unlisted please contact technical services.
Storage	Store at 4C in the dark.
Clonality	Monoclonal
Clone	SYP/7976R
Preservative	0.05% Sodium Azide
Isotype	IgG Kappa
Conjugate	DyLight 650
Purity	Protein A or G purified
Buffer	50mM Sodium Borate
Product Description	
Host	Rabbit
Gene ID	6855
Gene Symbol	SYP
Species	Human
Specificity/Sensitivity	This recombinant rabbit monoclonal antibody recognizes a protein of 38kDa that is identified as synaptophysin It is an N-glycosylated integral membrane protein found in neurons and endocrine cells Synaptophysin contains four transmembrane domains and may function as a gap junction-like channel This antibody identifies normal neuroendocrine cells and neuroendocrine neoplasms Diffuse, finely granular, cytoplasmic staining is observed, which probably correlates with the distribution of the antigen within neurosecretory vesicles Synaptophysin is an independent, broad-range marker of neural and neuroendocrine differentiation
Immunogen	Recombinant fragment (around aa 274-313) of human Synaptophysin protein (exact sequence is proprietary)
Notes	DyLight (R) is a trademark of Thermo Fisher Scientific Inc. and its subsidiaries.
Product Application Details	
Applications	Immunohistochemistry-Paraffin
Recommended Dilutions	Immunohistochemistry-Paraffin
Application Notes	Optimal dilution of this antibody should be experimentally determined.



Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA

Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
novus@novusbio.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada

Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom

Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: technical@novusbio.com
Orders: orders@novusbio.com
General: novus@novusbio.com

Products Related to NBP3-20998C

H00006855-P01-10ug	Recombinant Human Synaptophysin GST (N-Term) Protein
DBD00	BDNF [HRP]
NBP2-25170PEP	Synaptophysin Antibody Blocking Peptide
NB300-213	MAP2 Antibody

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/NBP3-20998C

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

