

Product Datasheet

Synaptophysin Antibody (SYP/7976R) [DyLight 550] NBP3-20998R

Unit Size: 0.1 ml

Store at 4C in the dark.

www.novusbio.com



technical@novusbio.com

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/NBP3-20998R

Updated 10/26/2023 v.20.1

Earn rewards for product
reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/NBP3-20998R



NBP3-20998R

Synaptophysin Antibody (SYP/7976R) [DyLight 550]

| Product Information | |
|------------------------------------|---|
| Unit Size | 0.1 ml |
| Concentration | Please see the vial label for concentration. If unlisted please contact technical services. |
| Storage | Store at 4C in the dark. |
| Clonality | Monoclonal |
| Clone | SYP/7976R |
| Preservative | 0.05% Sodium Azide |
| Isotype | IgG Kappa |
| Conjugate | DyLight 550 |
| Purity | Protein A or G purified |
| Buffer | 50mM Sodium Borate |
| Product Description | |
| Host | Rabbit |
| Gene ID | 6855 |
| Gene Symbol | SYP |
| Species | Human |
| Specificity/Sensitivity | This recombinant rabbit monoclonal antibody recognizes a protein of 38kDa that is identified as synaptophysin It is an N-glycosylated integral membrane protein found in neurons and endocrine cells Synaptophysin contains four transmembrane domains and may function as a gap junction-like channel This antibody identifies normal neuroendocrine cells and neuroendocrine neoplasms Diffuse, finely granular, cytoplasmic staining is observed, which probably correlates with the distribution of the antigen within neurosecretory vesicles Synaptophysin is an independent, broad-range marker of neural and neuroendocrine differentiation |
| Immunogen | Recombinant fragment (around aa 274-313) of human Synaptophysin protein (exact sequence is proprietary) |
| Notes | DyLight (R) is a trademark of Thermo Fisher Scientific Inc. and its subsidiaries. |
| Product Application Details | |
| Applications | Immunohistochemistry-Paraffin |
| Recommended Dilutions | Immunohistochemistry-Paraffin |
| Application Notes | Optimal dilution of this antibody should be experimentally determined. |



Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA

Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
novus@novusbio.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada

Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom

Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: technical@novusbio.com
Orders: orders@novusbio.com
General: novus@novusbio.com

Products Related to NBP3-20998R

| | |
|--------------------|--|
| H00006855-P01-10ug | Recombinant Human Synaptophysin GST (N-Term) Protein |
| DBD00 | BDNF [HRP] |
| NBP2-25170PEP | Synaptophysin Antibody Blocking Peptide |
| NB300-213 | MAP2 Antibody |

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/NBP3-20998R

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

