

# Product Datasheet

## ADAL Antibody NBP3-21262-100ul

Unit Size: 100 ul

Store at 4C short term. Aliquot and store at -20C long term. Avoid freeze-thaw cycles.

[www.novusbio.com](http://www.novusbio.com)



[technical@novusbio.com](mailto:technical@novusbio.com)

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:  
[www.novusbio.com/NBP3-21262](http://www.novusbio.com/NBP3-21262)

Updated 5/3/2023 v.20.1

Earn rewards for product  
reviews and publications.

Submit a publication at [www.novusbio.com/publications](http://www.novusbio.com/publications)

Submit a review at [www.novusbio.com/reviews/destination/NBP3-21262](http://www.novusbio.com/reviews/destination/NBP3-21262)



**NBP3-21262-100ul**

ADAL Antibody

**Product Information**

|                      |                                                                                                                  |
|----------------------|------------------------------------------------------------------------------------------------------------------|
| <b>Unit Size</b>     | 100 ul                                                                                                           |
| <b>Concentration</b> | Concentrations vary lot to lot. See vial label for concentration. If unlisted please contact technical services. |
| <b>Storage</b>       | Store at 4C short term. Aliquot and store at -20C long term. Avoid freeze-thaw cycles.                           |
| <b>Clonality</b>     | Polyclonal                                                                                                       |
| <b>Preservative</b>  | 0.02% Sodium Azide                                                                                               |
| <b>Isotype</b>       | IgG                                                                                                              |
| <b>Purity</b>        | Affinity purified                                                                                                |
| <b>Buffer</b>        | PBS, pH 7.2, 40% glycerol                                                                                        |

**Product Description**


|                    |                                                                                                                                                                                  |
|--------------------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Host</b>        | Rabbit                                                                                                                                                                           |
| <b>Gene ID</b>     | 161823                                                                                                                                                                           |
| <b>Gene Symbol</b> | ADAL                                                                                                                                                                             |
| <b>Species</b>     | Human                                                                                                                                                                            |
| <b>Immunogen</b>   | This antibody was developed against Recombinant Protein corresponding to amino acids:<br>DILMVTKDVIKEFADDGVKYLELRSTPRRENATGMTKKTYVESILEGIKQSKQEN<br>LDIDVRYLIAVDRRGGPLVAKETVKLAE |

**Product Application Details**

|                              |                                                         |
|------------------------------|---------------------------------------------------------|
| <b>Applications</b>          | Immunocytochemistry/Immunofluorescence                  |
| <b>Recommended Dilutions</b> | Immunocytochemistry/Immunofluorescence 0.25-2 ug/ml     |
| <b>Application Notes</b>     | ICC/IF Fixation Permeabilization: Use PFA/Triton X-100. |

**Images**

Staining of human cell line RH-30 shows localization to nucleoplasm & vesicles.





### **Novus Biologicals USA**

10730 E. Briarwood Avenue  
Centennial, CO 80112  
USA  
Phone: 303.730.1950  
Toll Free: 1.888.506.6887  
Fax: 303.730.1966  
novus@novusbio.com

### **Bio-Techne Canada**

21 Canmotor Ave  
Toronto, ON M8Z 4E6  
Canada  
Phone: 905.827.6400  
Toll Free: 855.668.8722  
Fax: 905.827.6402  
canada.inquires@bio-techne.com

### **Bio-Techne Ltd**

19 Barton Lane  
Abingdon Science Park  
Abingdon, OX14 3NB, United Kingdom  
Phone: (44) (0) 1235 529449  
Free Phone: 0800 37 34 15  
Fax: (44) (0) 1235 533420  
info.EMEA@bio-techne.com

### **General Contact Information**

www.novusbio.com  
Technical Support: technical@novusbio.com  
Orders: orders@novusbio.com  
General: novus@novusbio.com

### **Products Related to NBP3-21262-100ul**

---

|                  |                                               |
|------------------|-----------------------------------------------|
| HAF008           | Goat anti-Rabbit IgG Secondary Antibody [HRP] |
| NB7156           | Goat anti-Rabbit IgG (H+L) Secondary Antibody |
| NBP2-24891       | Rabbit IgG Isotype Control                    |
| NBP2-51649-0.1mg | Recombinant Human ADAL His Protein            |

---

### **Limitations**

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit [www.novusbio.com/guarantee](http://www.novusbio.com/guarantee)

Earn gift cards/discounts by submitting a review: [www.novusbio.com/reviews/submit/NBP3-21262](http://www.novusbio.com/reviews/submit/NBP3-21262)

Earn gift cards/discounts by submitting a publication using this product:  
[www.novusbio.com/publications](http://www.novusbio.com/publications)

