

Product Datasheet

DGCR8 Recombinant Protein Antigen NBP3-21419PEP

Unit Size: 100 ul

Store at -20C. Avoid freeze-thaw cycles.

www.novusbio.com



technical@novusbio.com

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/NBP3-21419PEP

Updated 6/13/2023 v.20.1

Earn rewards for product
reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/NBP3-21419PEP



NBP3-21419PEP**DGCR8 Recombinant Protein Antigen**

Product Information	
Unit Size	100 ul
Concentration	Please see the vial label for concentration. If unlisted please contact technical services.
Storage	Store at -20C. Avoid freeze-thaw cycles.
Preservative	No Preservative
Purity	>80% by SDS-PAGE and Coomassie blue staining
Buffer	PBS and 1M Urea, pH 7.4.
Target Molecular Weight	37 kDa
Product Description	
Description	<p>A recombinant protein antigen with a N-terminal His6-ABP tag corresponding to human DGCR8</p> <p>Source: <i>E.coli</i></p> <p>Amino Acid Sequence: YGASLLSKGSFSKGRLLIDPNCSGHSPRTARHAPAVRKFSFDLKLLKDVKISVS FTESCRSKDRKVLTYGAERDVRAECGLLLSPVSGDVHACPFGGSVGDGVGIG GESADKKDEENELDQEKRVEYAVLDELEDFTDNLELDEEGAGGFTAIVQRD RVDEEALNFPYEDDFDNDVDAL</p> <p>Fusion Tag: N-terminal His6ABP (ABP = Albumin Binding Protein derived from Streptococcal Protein G)</p> <p>This product is intended to be used as a blocking antigen for antibody competition assays. Any other use of this antigen is done at the risk of the user. The use of this product for commercial production is strictly prohibited. Please contact technical support if you have any questions.</p>
Gene ID	54487
Gene Symbol	DGCR8
Species	Human
Product Application Details	
Applications	Antibody Competition
Recommended Dilutions	Antibody Competition 10-100 molar excess
Application Notes	<p>This recombinant antigen is only intended to be used as a blocking agent to confirm antibody specificity with the corresponding antibody, catalog number NBP3-21419. It is purified by IMAC chromatography, and the expected concentration is greater than 0.5 mg/ml. For current lot information, including availability, please contact our technical support team click nb-technical@bio-technique.com. For further blocking peptide related information and a protocol, click here.</p>



Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA

Phone: 303.730.1950

Toll Free: 1.888.506.6887

Fax: 303.730.1966

nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada

Phone: 905.827.6400

Toll Free: 855.668.8722

Fax: 905.827.6402

canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom

Phone: (44) (0) 1235 529449

Free Phone: 0800 37 34 15

Fax: (44) (0) 1235 533420

info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com

Technical Support: nb-technical@bio-techne.com

Orders: nb-customerservice@bio-techne.com

General: novus@novusbio.com

Products Related to NBP3-21419PEP

H00054487-Q01-10ug	Recombinant Human DGCR8 GST (N-Term) Protein
NBP1-28751	GW182 Antibody
H00054487-M01	DGCR8 Antibody (4G3)
NBP1-71691-0.1mg	Dicer Antibody - BSA Free

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Peptides and proteins are guaranteed for 3 months from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/NBP3-21419PEP

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

