

# Product Datasheet

## RAR gamma/NR1B3 Antibody (SR1891) - BSA Free NBP3-21974-100ul

Unit Size: 100 ul

Store at 4C short term. Aliquot and store at -20C long term. Avoid freeze-thaw cycles.

[www.novusbio.com](http://www.novusbio.com)



[technical@novusbio.com](mailto:technical@novusbio.com)

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:  
[www.novusbio.com/NBP3-21974](http://www.novusbio.com/NBP3-21974)

Updated 9/9/2025 v.20.1

Earn rewards for product  
reviews and publications.

Submit a publication at [www.novusbio.com/publications](http://www.novusbio.com/publications)

Submit a review at [www.novusbio.com/reviews/destination/NBP3-21974](http://www.novusbio.com/reviews/destination/NBP3-21974)



**NBP3-21974-100ul**

RAR gamma/NR1B3 Antibody (SR1891) - BSA Free


| Product Information     |                                                                                             |
|-------------------------|---------------------------------------------------------------------------------------------|
| Unit Size               | 100 ul                                                                                      |
| Concentration           | Please see the vial label for concentration. If unlisted please contact technical services. |
| Storage                 | Store at 4C short term. Aliquot and store at -20C long term. Avoid freeze-thaw cycles.      |
| Clonality               | Monoclonal                                                                                  |
| Clone                   | SR1891                                                                                      |
| Preservative            | 0.02% Sodium Azide                                                                          |
| Isotype                 | IgG                                                                                         |
| Purity                  | Affinity purified                                                                           |
| Buffer                  | PBS, pH 7.4, 150mM NaCl, 50% glycerol.                                                      |
| Target Molecular Weight | 50 kDa                                                                                      |

| Product Description     |                                                                                                                                                                                                                                 |
|-------------------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Description             | Novus Biologicals Rabbit RAR gamma/NR1B3 Antibody (SR1891) - BSA Free (NBP3-21974) is a recombinant monoclonal antibody validated for use in WB and ICC/IF. All Novus Biologicals antibodies are covered by our 100% guarantee. |
| Host                    | Rabbit                                                                                                                                                                                                                          |
| Gene ID                 | 5916                                                                                                                                                                                                                            |
| Gene Symbol             | RARG                                                                                                                                                                                                                            |
| Species                 | Human, Rat                                                                                                                                                                                                                      |
| Specificity/Sensitivity | Detects endogenous levels of total RAR gamma/NR1B3                                                                                                                                                                              |
| Immunogen               | A synthesized peptide derived from human RAR gamma/NR1B3 (Uniprot #: P13631)                                                                                                                                                    |

| Product Application Details |                                                                               |
|-----------------------------|-------------------------------------------------------------------------------|
| Applications                | Western Blot, Immunocytochemistry/ Immunofluorescence                         |
| Recommended Dilutions       | Western Blot 1:500-1:2000, Immunocytochemistry/ Immunofluorescence 1:50-1:200 |

**Images**

Western Blot: RAR gamma/NR1B3 Antibody (SR1891) [NBP3-21974] - Western blot analysis using NBP3-21974 in HeLa cell lysate.





### **Novus Biologicals USA**

10730 E. Briarwood Avenue  
Centennial, CO 80112  
USA  
Phone: 303.730.1950  
Toll Free: 1.888.506.6887  
Fax: 303.730.1966  
nb-customerservice@bio-techne.com

### **Bio-Techne Canada**

21 Canmotor Ave  
Toronto, ON M8Z 4E6  
Canada  
Phone: 905.827.6400  
Toll Free: 855.668.8722  
Fax: 905.827.6402  
canada.inquires@bio-techne.com

### **Bio-Techne Ltd**

19 Barton Lane  
Abingdon Science Park  
Abingdon, OX14 3NB, United Kingdom  
Phone: (44) (0) 1235 529449  
Free Phone: 0800 37 34 15  
Fax: (44) (0) 1235 533420  
info.EMEA@bio-techne.com

### **General Contact Information**

[www.novusbio.com](http://www.novusbio.com)  
Technical Support: nb-technical@bio-techne.com  
Orders: nb-customerservice@bio-techne.com  
General: novus@novusbio.com

### **Products Related to NBP3-21974-100ul**

|                    |                                                        |
|--------------------|--------------------------------------------------------|
| HAF008             | Goat anti-Rabbit IgG Secondary Antibody [HRP]          |
| NB7160             | Goat anti-Rabbit IgG (H+L) Secondary Antibody [HRP]    |
| NBP2-24891         | Rabbit IgG Isotype Control                             |
| H00005916-Q01-10ug | Recombinant Human RAR gamma/NR1B3 GST (N-Term) Protein |

### **Limitations**

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit [www.novusbio.com/guarantee](http://www.novusbio.com/guarantee)

Earn gift cards/discounts by submitting a review: [www.novusbio.com/reviews/submit/NBP3-21974](http://www.novusbio.com/reviews/submit/NBP3-21974)

Earn gift cards/discounts by submitting a publication using this product:  
[www.novusbio.com/publications](http://www.novusbio.com/publications)

