

# Product Datasheet

## PDE2A Antibody (SR1040)

### NBP3-22378-100ul

Unit Size: 100 ul

Store at 4C short term. Aliquot and store at -20C long term. Avoid freeze-thaw cycles.

[www.novusbio.com](http://www.novusbio.com)



[technical@novusbio.com](mailto:technical@novusbio.com)


Protocols, Publications, Related Products, Reviews, Research Tools and Images at:  
[www.novusbio.com/NBP3-22378](http://www.novusbio.com/NBP3-22378)

Updated 4/15/2024 v.20.1

Earn rewards for product  
reviews and publications.

Submit a publication at [www.novusbio.com/publications](http://www.novusbio.com/publications)

Submit a review at [www.novusbio.com/reviews/destination/NBP3-22378](http://www.novusbio.com/reviews/destination/NBP3-22378)



**NBP3-22378-100ul**

PDE2A Antibody (SR1040)


| Product Information     |                                                                                             |
|-------------------------|---------------------------------------------------------------------------------------------|
| Unit Size               | 100 ul                                                                                      |
| Concentration           | Please see the vial label for concentration. If unlisted please contact technical services. |
| Storage                 | Store at 4C short term. Aliquot and store at -20C long term. Avoid freeze-thaw cycles.      |
| Clonality               | Monoclonal                                                                                  |
| Clone                   | SR1040                                                                                      |
| Preservative            | 0.02% Sodium Azide                                                                          |
| Isotype                 | IgG                                                                                         |
| Purity                  | Affinity purified                                                                           |
| Buffer                  | PBS, pH 7.4, 150mM NaCl, 50% glycerol.                                                      |
| Target Molecular Weight | 106 kDa                                                                                     |

| Product Description     |                                                                    |
|-------------------------|--------------------------------------------------------------------|
| Host                    | Rabbit                                                             |
| Gene Symbol             | PDE2A                                                              |
| Species                 | Human, Mouse, Rat                                                  |
| Specificity/Sensitivity | Detects endogenous levels of total PDE2A                           |
| Immunogen               | A synthesized peptide derived from human PDE2A (Uniprot #: O00408) |

| Product Application Details |                           |
|-----------------------------|---------------------------|
| Applications                | Western Blot              |
| Recommended Dilutions       | Western Blot 1:500-1:2000 |

**Images**

Western Blot: PDE2A Antibody (SR1040) [NBP3-22378] - Western blot analysis using NBP3-22378 in Human fetal brain lysate.





### **Novus Biologicals USA**

10730 E. Briarwood Avenue  
Centennial, CO 80112  
USA  
Phone: 303.730.1950  
Toll Free: 1.888.506.6887  
Fax: 303.730.1966  
nb-customerservice@bio-techne.com

### **Bio-Techne Canada**

21 Canmotor Ave  
Toronto, ON M8Z 4E6  
Canada  
Phone: 905.827.6400  
Toll Free: 855.668.8722  
Fax: 905.827.6402  
canada.inquires@bio-techne.com

### **Bio-Techne Ltd**

19 Barton Lane  
Abingdon Science Park  
Abingdon, OX14 3NB, United Kingdom  
Phone: (44) (0) 1235 529449  
Free Phone: 0800 37 34 15  
Fax: (44) (0) 1235 533420  
info.EMEA@bio-techne.com

### **General Contact Information**

[www.novusbio.com](http://www.novusbio.com)  
Technical Support: nb-technical@bio-techne.com  
Orders: nb-customerservice@bio-techne.com  
General: novus@novusbio.com

### **Products Related to NBP3-22378-100ul**

|                    |                                                     |
|--------------------|-----------------------------------------------------|
| HAF008             | Goat anti-Rabbit IgG Secondary Antibody [HRP]       |
| NB7160             | Goat anti-Rabbit IgG (H+L) Secondary Antibody [HRP] |
| NBP2-24891         | Rabbit IgG Isotype Control                          |
| H00005138-Q01-10ug | Recombinant Human PDE2A GST (N-Term) Protein        |

### **Limitations**

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit [www.novusbio.com/guarantee](http://www.novusbio.com/guarantee)

Earn gift cards/discounts by submitting a review: [www.novusbio.com/reviews/submit/NBP3-22378](http://www.novusbio.com/reviews/submit/NBP3-22378)

Earn gift cards/discounts by submitting a publication using this product:  
[www.novusbio.com/publications](http://www.novusbio.com/publications)

