

Product Datasheet

beta-Actin Antibody NBP3-23319

Unit Size: 0.1 ml

Store at -20C.

www.novusbio.com



technical@novusbio.com

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/NBP3-23319

Updated 9/11/2024 v.20.1

Earn rewards for product
reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/NBP3-23319



NBP3-23319


beta-Actin Antibody

Product Information	
Unit Size	0.1 ml
Concentration	Please see the vial label for concentration. If unlisted please contact technical services.
Storage	Store at -20C.
Clonality	Polyclonal
Preservative	0.05% Sodium Azide
Isotype	IgG
Purity	Antigen Affinity-purified
Buffer	PBS, 1 mg/ml BSA and 50% glycerol
Target Molecular Weight	42 kDa
Product Description	
Description	Stable for at least 1 year at -20°C.
Host	Rabbit
Gene ID	60
Gene Symbol	ACTB
Species	Human, Mouse, Rat, Chicken
Specificity/Sensitivity	This antibody detects a 42 kDa protein corresponding to the molecular mass of Actin on SDS-PAGE immunoblots of human A431, SYF, and HUVEC cells, as well as mouse C2C12 cells.
Immunogen	Actin synthetic peptide (coupled to KLH) corresponding to amino acid residues in the N-terminal region of human beta-Actin. This sequence is identical to similar regions in all four alpha-actins, as well as in gamma-actin, and is well conserved in actins from most eukaryotic species. [UniProt# P60709]
Product Application Details	
Applications	Western Blot, Immunocytochemistry/ Immunofluorescence, Immunohistochemistry, Immunoprecipitation
Recommended Dilutions	Western Blot 1:1000, Immunohistochemistry 1:100, Immunocytochemistry/ Immunofluorescence 1:50, Immunoprecipitation 1:25




Images

Immunocytochemistry/Immunofluorescence: beta-Actin Antibody [NBP3-23319] - Immunocytochemical labeling using anti-Actin (N-terminal) and anti-Actin (Tyr-53) polyclonal antibodies in C2C12 cells control (left) or treated with pervanadate (1 mM) for 30 min (middle). The cells were fixed in paraformaldehyde and permeabilized in acetone. Both antibodies were used in the presence of blocking peptide: Actin (N-terminal) peptide or phospho-Actin (Tyr-53) peptide, respectively (right).



Western Blot: beta-Actin Antibody [NBP3-23319] - Western blot analysis of mouse C2C12 cells untreated (lanes 1 & 3), or treated with pervanadate (1 mM) for 30 min (lanes 2 & 4). The blot was probed with anti-Actin (N-terminal) antibody (lanes 1 & 2) or anti-Actin (Tyr-53) antibody (lanes 3 & 4).





Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: nb-technical@bio-techne.com
Orders: nb-customerservice@bio-techne.com
General: novus@novusbio.com

Products Related to NBP3-23319

HAF008	Goat anti-Rabbit IgG Secondary Antibody [HRP]
NB7160	Goat anti-Rabbit IgG (H+L) Secondary Antibody [HRP]
NBP2-24891	Rabbit IgG Isotype Control
NB600-503PEP	beta-Actin Antibody Blocking Peptide

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/NBP3-23319

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

