

Product Datasheet

beta-Actin Antibody (M008)

NBP3-23324

Unit Size: 0.1 ml

Store at -20C.

www.novusbio.com



technical@novusbio.com

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/NBP3-23324

Updated 9/11/2024 v.20.1

Earn rewards for product
reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/NBP3-23324



NBP3-23324


beta-Actin Antibody (M008)

Product Information	
Unit Size	0.1 ml
Concentration	Please see the vial label for concentration. If unlisted please contact technical services.
Storage	Store at -20C.
Clonality	Monoclonal
Clone	M008
Preservative	0.05% Sodium Azide
Isotype	IgG2b
Purity	Protein G purified
Buffer	PBS, 1 mg/ml BSA and 50% glycerol
Target Molecular Weight	42 kDa
Product Description	
Description	Stable for at least 1 year at -20°C.
Host	Mouse
Gene ID	60
Gene Symbol	ACTB
Species	Human
Specificity/Sensitivity	Clone M008 detects a 42 kDa* protein corresponding to the molecular mass of beta-actin on SDS-PAGE immunoblots of human cancer cell lines, as well as human recombinant beta actin. This actin antibody preferentially detects human beta-actin with only weak reactivity toward actins in rat, mouse, or rabbit.
Immunogen	Clone (M008) was generated from a proprietary antigen related to human beta-Actin in MDA-MB-231 breast cancer cell line. [UniProt# P60709]
Product Application Details	
Applications	Western Blot, Immunocytochemistry/ Immunofluorescence, Immunoprecipitation
Recommended Dilutions	Western Blot 1:1000, Immunocytochemistry/ Immunofluorescence 1:50, Immunoprecipitation 1:100




Images

Western Blot: beta-Actin Antibody (M008) [NBP3-23324] - Western blot analysis of human HUVEC-CS (lane 1), rabbit spleen fibroblast (lane 2), human Jurkat (lane 3), human LNCaP (lane 4), human HeLa (lane 5), and mouse F9 (lane 6) cell lysates. The blot was probed with mouse monoclonal anti-beta-Actin (NBP3-23324) at 1:1000 (lanes 1-6).



Immunocytochemistry/Immunofluorescence: beta-Actin Antibody (M008) [NBP3-23324] - Immunocytochemical labeling of beta-Actin in paraformaldehyde fixed human MeWo cells. The cells were labeled with mouse monoclonal anti-beta-Actin (clone M008). The antibody was detected using goat anti-mouse DyLight (R) 594.





Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: nb-technical@bio-techne.com
Orders: nb-customerservice@bio-techne.com
General: novus@novusbio.com

Products Related to NBP3-23324

HAF007	Goat anti-Mouse IgG Secondary Antibody [HRP]
NB720-B	Rabbit anti-Mouse IgG (H+L) Secondary Antibody [Biotin]
NBP2-27231	Mouse IgG2b Isotype Control (MPC-11)
NB600-503PEP	beta-Actin Antibody Blocking Peptide

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/NBP3-23324

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

