

# Product Datasheet

## **FHOD1 Antibody (M352) NBP3-23340**

Unit Size: 0.1 ml

Store at -20C.

[www.novusbio.com](http://www.novusbio.com)



[technical@novusbio.com](mailto:technical@novusbio.com)

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:  
[www.novusbio.com/NBP3-23340](http://www.novusbio.com/NBP3-23340)

Updated 8/28/2023 v.20.1

**Earn rewards for product  
reviews and publications.**

Submit a publication at [www.novusbio.com/publications](http://www.novusbio.com/publications)

Submit a review at [www.novusbio.com/reviews/destination/NBP3-23340](http://www.novusbio.com/reviews/destination/NBP3-23340)



**NBP3-23340****FHOD1 Antibody (M352)**

| <b>Product Information</b>         |                                                                                                                                                                                                 |
|------------------------------------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Unit Size</b>                   | 0.1 ml                                                                                                                                                                                          |
| <b>Concentration</b>               | Please see the vial label for concentration. If unlisted please contact technical services.                                                                                                     |
| <b>Storage</b>                     | Store at -20C.                                                                                                                                                                                  |
| <b>Clonality</b>                   | Monoclonal                                                                                                                                                                                      |
| <b>Clone</b>                       | M352                                                                                                                                                                                            |
| <b>Preservative</b>                | 0.05% Sodium Azide                                                                                                                                                                              |
| <b>Isotype</b>                     | IgG1                                                                                                                                                                                            |
| <b>Purity</b>                      | Protein A purified                                                                                                                                                                              |
| <b>Buffer</b>                      | PBS, 1 mg/ml BSA and 50% glycerol                                                                                                                                                               |
| <b>Target Molecular Weight</b>     | 140 kDa                                                                                                                                                                                         |
| <b>Product Description</b>         |                                                                                                                                                                                                 |
| <b>Description</b>                 | Stable for at least 1 year at -20°C.                                                                                                                                                            |
| <b>Host</b>                        | Mouse                                                                                                                                                                                           |
| <b>Gene ID</b>                     | 29109                                                                                                                                                                                           |
| <b>Gene Symbol</b>                 | FHOD1                                                                                                                                                                                           |
| <b>Species</b>                     | Human                                                                                                                                                                                           |
| <b>Specificity/Sensitivity</b>     | The antibody detects a 140 kDa protein corresponding to the apparent molecular mass of FHOD1 on SDS-PAGE immunoblots of human K562 and Jurkat, and has weak reactivity to FHOD1 in mouse C2C12. |
| <b>Immunogen</b>                   | FHOD1 recombinant protein that included the full length human FHOD1 sequence. This sequence has 50% homology to human FHOD3, and low homology to other formin family members. [UniProt# Q9Y613] |
| <b>Product Application Details</b> |                                                                                                                                                                                                 |
| <b>Applications</b>                | Western Blot, Immunohistochemistry                                                                                                                                                              |
| <b>Recommended Dilutions</b>       | Western Blot 1:250, Immunohistochemistry 1:50                                                                                                                                                   |





### **Novus Biologicals USA**

10730 E. Briarwood Avenue  
Centennial, CO 80112  
USA

Phone: 303.730.1950  
Toll Free: 1.888.506.6887  
Fax: 303.730.1966  
novus@novusbio.com

### **Bio-Techne Canada**

21 Canmotor Ave  
Toronto, ON M8Z 4E6  
Canada  
Phone: 905.827.6400  
Toll Free: 855.668.8722  
Fax: 905.827.6402  
canada.inquires@bio-techne.com

### **Bio-Techne Ltd**

19 Barton Lane  
Abingdon Science Park  
Abingdon, OX14 3NB, United Kingdom  
Phone: (44) (0) 1235 529449  
Free Phone: 0800 37 34 15  
Fax: (44) (0) 1235 533420  
info.EMEA@bio-techne.com

### **General Contact Information**

www.novusbio.com  
Technical Support: technical@novusbio.com  
Orders: orders@novusbio.com  
General: novus@novusbio.com

### **Products Related to NBP3-23340**

|                    |                                                         |
|--------------------|---------------------------------------------------------|
| HAF007             | Goat anti-Mouse IgG Secondary Antibody [HRP]            |
| NB720-B            | Rabbit anti-Mouse IgG (H+L) Secondary Antibody [Biotin] |
| NBP1-97005-0.5mg   | Mouse IgG1 Isotype Control (MG1)                        |
| H00029109-Q01-10ug | Recombinant Human FHOD1 GST (N-Term) Protein            |

### **Limitations**

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit [www.novusbio.com/guarantee](http://www.novusbio.com/guarantee)

Earn gift cards/discounts by submitting a review: [www.novusbio.com/reviews/submit/NBP3-23340](http://www.novusbio.com/reviews/submit/NBP3-23340)

Earn gift cards/discounts by submitting a publication using this product:  
[www.novusbio.com/publications](http://www.novusbio.com/publications)

