

# Product Datasheet

## Involucrin Antibody (rIVRN/9323) [DyLight 350] NBP3-24321UV

Unit Size: 0.1 ml

Store at 4C in the dark.

[www.novusbio.com](http://www.novusbio.com)



[technical@novusbio.com](mailto:technical@novusbio.com)

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:  
[www.novusbio.com/NBP3-24321UV](http://www.novusbio.com/NBP3-24321UV)

Updated 11/18/2023 v.20.1

Earn rewards for product  
reviews and publications.

Submit a publication at [www.novusbio.com/publications](http://www.novusbio.com/publications)

Submit a review at [www.novusbio.com/reviews/destination/NBP3-24321UV](http://www.novusbio.com/reviews/destination/NBP3-24321UV)



**NBP3-24321UV**

Involucrin Antibody (rIVRN/9323) [DyLight 350]

| Product Information         |                                                                                                                                                                                                                                                                                                                                                                                                                            |
|-----------------------------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Unit Size                   | 0.1 ml                                                                                                                                                                                                                                                                                                                                                                                                                     |
| Concentration               | Please see the vial label for concentration. If unlisted please contact technical services.                                                                                                                                                                                                                                                                                                                                |
| Storage                     | Store at 4C in the dark.                                                                                                                                                                                                                                                                                                                                                                                                   |
| Clonality                   | Monoclonal                                                                                                                                                                                                                                                                                                                                                                                                                 |
| Clone                       | rIVRN/9323                                                                                                                                                                                                                                                                                                                                                                                                                 |
| Preservative                | 0.05% Sodium Azide                                                                                                                                                                                                                                                                                                                                                                                                         |
| Isotype                     | IgG2b Kappa                                                                                                                                                                                                                                                                                                                                                                                                                |
| Conjugate                   | DyLight 350                                                                                                                                                                                                                                                                                                                                                                                                                |
| Purity                      | Protein A purified                                                                                                                                                                                                                                                                                                                                                                                                         |
| Buffer                      | 50mM Sodium Borate                                                                                                                                                                                                                                                                                                                                                                                                         |
| Product Description         |                                                                                                                                                                                                                                                                                                                                                                                                                            |
| Host                        | Mouse                                                                                                                                                                                                                                                                                                                                                                                                                      |
| Gene ID                     | 3713                                                                                                                                                                                                                                                                                                                                                                                                                       |
| Gene Symbol                 | IVL                                                                                                                                                                                                                                                                                                                                                                                                                        |
| Species                     | Human                                                                                                                                                                                                                                                                                                                                                                                                                      |
| Specificity/Sensitivity     | It recognizes a protein of 66kDa-170kDa, identified as involucrin. In Western blotting of cultured human keratinocytes, this monoclonal antibody reacts with a 120kDa protein. It stains the involucrin in a variety of sizes: 170kDa in MCF-7 cells, a doublet of ~115kDa and 150kDa in gorilla and owl monkey, 66kDa in dog, and a doublet of 105kDa in pig. Its epitope maps between codon 421-568 of human involucrin. |
| Immunogen                   | Recombinant human full-length Involucrin protein (IVRN)                                                                                                                                                                                                                                                                                                                                                                    |
| Notes                       | DyLight (R) is a trademark of Thermo Fisher Scientific Inc. and its subsidiaries.                                                                                                                                                                                                                                                                                                                                          |
| Product Application Details |                                                                                                                                                                                                                                                                                                                                                                                                                            |
| Applications                | Immunohistochemistry-Paraffin                                                                                                                                                                                                                                                                                                                                                                                              |
| Recommended Dilutions       | Immunohistochemistry-Paraffin                                                                                                                                                                                                                                                                                                                                                                                              |
| Application Notes           | Optimal dilution of this antibody should be experimentally determined.                                                                                                                                                                                                                                                                                                                                                     |





### **Novus Biologicals USA**

10730 E. Briarwood Avenue  
Centennial, CO 80112  
USA  
Phone: 303.730.1950  
Toll Free: 1.888.506.6887  
Fax: 303.730.1966  
nb-customerservice@bio-techne.com

### **Bio-Techne Canada**

21 Canmotor Ave  
Toronto, ON M8Z 4E6  
Canada  
Phone: 905.827.6400  
Toll Free: 855.668.8722  
Fax: 905.827.6402  
canada.inquires@bio-techne.com

### **Bio-Techne Ltd**

19 Barton Lane  
Abingdon Science Park  
Abingdon, OX14 3NB, United Kingdom  
Phone: (44) (0) 1235 529449  
Free Phone: 0800 37 34 15  
Fax: (44) (0) 1235 533420  
info.EMEA@bio-techne.com

### **General Contact Information**

[www.novusbio.com](http://www.novusbio.com)  
Technical Support: nb-technical@bio-techne.com  
Orders: nb-customerservice@bio-techne.com  
General: novus@novusbio.com

### **Products Related to NBP3-24321UV**

---

|               |                                                                    |
|---------------|--------------------------------------------------------------------|
| NBP1-43317UV  | Mouse IgG2b Kappa Light Chain Isotype Control (MG2b) [DyLight 350] |
| NBP2-33742PEP | Involucrin Recombinant Protein Antigen                             |
| 236-EG-200    | EGF [Unconjugated]                                                 |
| NB100-355     | RPE65 Antibody (401.8B11.3D9) - BSA Free                           |

---

### **Limitations**

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit [www.novusbio.com/guarantee](http://www.novusbio.com/guarantee)

Earn gift cards/discounts by submitting a review: [www.novusbio.com/reviews/submit/NBP3-24321UV](http://www.novusbio.com/reviews/submit/NBP3-24321UV)

Earn gift cards/discounts by submitting a publication using this product:  
[www.novusbio.com/publications](http://www.novusbio.com/publications)

