

Product Datasheet

Human DEFA6 ELISA Kit (Colorimetric) NBP3-24605

Unit Size: 1 Kit

Storage of components varies. See protocol for specific instructions.

www.novusbio.com



technical@novusbio.com

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/NBP3-24605

Updated 7/30/2025 v.20.1

Earn rewards for product
reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/NBP3-24605



NBP3-24605**Human DEFA6 ELISA Kit (Colorimetric)**


| Product Information | |
|----------------------------|--|
| Unit Size | 1 Kit |
| Concentration | Concentration is not relevant for this product. Please see the protocols for proper use of this product. |
| Storage | Storage of components varies. See protocol for specific instructions. |

| Product Description | |
|--------------------------------|--|
| Description | Assay Length: 1-5 hours |
| Gene ID | 1671 |
| Gene Symbol | DEFA6 |
| Species | Human |
| Specificity/Sensitivity | This assay has high sensitivity and excellent specificity for detection of human DEFA6. No significant cross-reactivity or interference between human DEFA6 and analogues was observed. |
| Kit Components | Assay plate (12 x 8 coated Microwells), Standard (Freeze dried), Biotin-antibody (100 x concentrate), HRP-avidin (100 x concentrate), Biotin-antibody Diluent, HRP-avidin Diluent, Sample Diluent, Wash Buffer (25 x concentrate), TMB Substrate, Stop Solution, Adhesive Strip (For 96 wells), Instruction manual |
| Standard Curve Range | 0.156 ng/mL - 10 ng/mL |
| Sensitivity | 0.039 ng/mL |
| Inter Assay | %CV < 10.00 |
| Intra Assay | %CV < 8.00 |
| Assay Type | Sandwich ELISA |
| Spiking Recovery | 81 - 100% |
| Suitable Sample Type | Serum, plasma and tissue homogenates |
| Sample Volume | 50 - 100 μ L |

| Product Application Details | |
|------------------------------------|-------|
| Applications | ELISA |
| Recommended Dilutions | ELISA |

Images

Human DEFA6 ELISA Kit (Colorimetric) [NBP3-24605] - Standard Curve Reference





Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: nb-technical@bio-techne.com
Orders: nb-customerservice@bio-techne.com
General: novus@novusbio.com

Products Related to NBP3-24605

| | |
|--------------------|--|
| H00001671-P01-10ug | Recombinant Human DEFA6 GST (N-Term) Protein |
| MAB208-100 | CXCL8/IL-8 Antibody (6217) [Unconjugated] |
| NBP2-97543-100ul | DEFA6 Antibody - BSA Free |
| NB100-524 | NOD2 Antibody (2D9) - BSA Free |

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. ELISA Kits are guaranteed for 6 months from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/NBP3-24605

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

