

Product Datasheet

Galectin-3BP/MAC-2BP/LGALS3BP Antibody (HL2168) - Azide and BSA Free NBP3-25550

Unit Size: 100 ul

Store at 4C short term. Aliquot and store at -20C long term. Avoid freeze-thaw cycles.

www.novusbio.com



technical@novusbio.com

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/NBP3-25550

Updated 12/19/2023 v.20.1

Earn rewards for product
reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/NBP3-25550



NBP3-25550

Galectin-3BP/MAC-2BP/LGALS3BP Antibody (HL2168) - Azide and BSA Free

Product Information

Unit Size	100 ul
Concentration	Concentrations vary lot to lot. See vial label for concentration. If unlisted please contact technical services.
Storage	Store at 4C short term. Aliquot and store at -20C long term. Avoid freeze-thaw cycles.
Clonality	Monoclonal
Clone	HL2168
Preservative	No Preservative
Isotype	IgG
Purity	Protein A purified
Buffer	PBS

Product Description


Host	Rabbit
Gene ID	3959
Gene Symbol	LGALS3BP
Species	Human, Mouse
Immunogen	Recombinant fragment of human Galectin-3BP/MAC-2BP/LGALS3BP.

Product Application Details


Applications	Western Blot, Immunocytochemistry/ Immunofluorescence, Immunohistochemistry, Immunohistochemistry-Paraffin
Recommended Dilutions	Western Blot 1:500-1:3000, Immunohistochemistry, Immunocytochemistry/ Immunofluorescence, Immunohistochemistry-Paraffin

Images


Immunocytochemistry/Immunofluorescence: Galectin-3BP/MAC-2BP/LGALS3BP Antibody (HL2168) - Azide and BSA Free [NBP3-25550] - LGALS3BP antibody [HL2168] detects LGALS3BP protein at cytoplasmic vesicle by immunofluorescent analysis. Sample: A431 cells were fixed in ice-cold MeOH for 5 min. Green: LGALS3BP stained by LGALS3BP antibody [HL2168] (NBP3-25550) diluted at 1:500. Red: alpha Tubulin, a cytoskeleton marker, stained by alpha Tubulin antibody [GT114] (NBP2-43837) diluted at 1:1000. Blue: Fluoroshield with DAPI .




Immunohistochemistry-Paraffin: Galectin-3BP/MAC-2BP/LGALS3BP Antibody (HL2168) - Azide and BSA Free [NBP3-25550] - LGALS3BP antibody [HL2168] detects secreted LGALS3BP protein by immunohistochemical analysis. Sample: Paraffin-embedded mouse adrenal gland. LGALS3BP stained by LGALS3BP antibody [HL2168] (NBP3-25550) diluted at 1:100. Antigen Retrieval: Citrate buffer, pH 6.0, 15 min




Immunocytochemistry/Immunofluorescence: Galectin-3BP/MAC-2BP/LGALS3BP Antibody (HL2168) - Azide and BSA Free [NBP3-25550] - LGALS3BP antibody [HL2168] detects LGALS3BP protein at endoplasmic reticulum and cell junction by immunofluorescent analysis. Sample: HeLa cells were fixed in 4% paraformaldehyde at RT for 15 min. Green: LGALS3BP stained by LGALS3BP antibody [HL2168] (NBP3-25550) diluted at 1:500. Red: alpha Tubulin, a cytoskeleton marker, stained by alpha Tubulin antibody [GT114] (NBP2-43837) diluted at 1:1000. Blue: Fluoroshield with DAPI .



Western Blot: Galectin-3BP/MAC-2BP/LGALS3BP Antibody (HL2168) - Azide and BSA Free [NBP3-25550] - Various whole cell extracts (30 ug) were separated by 7.5% SDS-PAGE, and the membrane was blotted with LGALS3BP antibody [HL2168] (NBP3-25550) diluted at 1:1000. The HRP-conjugated anti-rabbit IgG antibody was used to detect the primary antibody. Corresponding RNA expression data for the same cell lines are based on Human Protein Atlas program.



Western Blot: Galectin-3BP/MAC-2BP/LGALS3BP Antibody (HL2168) - Azide and BSA Free [NBP3-25550] - Untreated (-) and treated (+) A431 whole cell extracts (30 ug) were separated by 7.5% SDS-PAGE, and the membrane was blotted with LGALS3BP antibody [HL2168] (NBP3-25550) diluted at 1:1000. The HRP-conjugated anti-rabbit IgG antibody was used to detect the primary antibody.





Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: nb-technical@bio-techne.com
Orders: nb-customerservice@bio-techne.com
General: novus@novusbio.com

Products Related to NBP3-25550

HAF008	Goat anti-Rabbit IgG Secondary Antibody [HRP]
NB7160	Goat anti-Rabbit IgG (H+L) Secondary Antibody [HRP]
NBP2-24891	Rabbit IgG Isotype Control
NBP1-89346PEP	Galectin-3BP/MAC-2BP/LGALS3BP Recombinant Protein Antigen

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/NBP3-25550

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

