

Product Datasheet

Respiratory Syncytial Virus Glycoprotein G Antibody (HL1733) - Azide and BSA Free NBP3-25633

Unit Size: 100 ul

Store at 4C short term. Aliquot and store at -20C long term. Avoid freeze-thaw cycles.

www.novusbio.com



technical@novusbio.com

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/NBP3-25633

Updated 9/9/2025 v.20.1

Earn rewards for product
reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/NBP3-25633



NBP3-25633**Respiratory Syncytial Virus Glycoprotein G Antibody (HL1733) - Azide and BSA Free**

Product Information	
Unit Size	100 ul
Concentration	Concentrations vary lot to lot. See vial label for concentration. If unlisted please contact technical services.
Storage	Store at 4C short term. Aliquot and store at -20C long term. Avoid freeze-thaw cycles.
Clonality	Monoclonal
Clone	HL1733
Preservative	No Preservative
Isotype	IgG
Purity	Protein A purified
Buffer	PBS


Product Description	
Description	Novus Biologicals Rabbit Respiratory Syncytial Virus Glycoprotein G Antibody (HL1733) - Azide and BSA Free (NBP3-25633) is a recombinant monoclonal antibody validated for use in IHC, WB and ICC/IF. All Novus Biologicals antibodies are covered by our 100% guarantee.
Host	Rabbit
Species	Virus
Reactivity Notes	Respiratory Syncytial Virus type B
Immunogen	Recombinant fragment of Respiratory Syncytial Virus Glycoprotein G protein.

Product Application Details	
Applications	Western Blot, Immunohistochemistry-Paraffin, Immunocytochemistry/Immunofluorescence, Immunohistochemistry
Recommended Dilutions	Western Blot 1:1000-1:10000, Immunohistochemistry, Immunocytochemistry/Immunofluorescence, Immunohistochemistry-Paraffin




Images


Western Blot: Respiratory Syncytial Virus Glycoprotein G Antibody (HL1733) - Azide and BSA Free [NBP3-25633] - Non-transfected (-) and transfected (+) 293T whole cell extracts (30 ug) were separated by 12% SDS-PAGE, and the membrane was blotted with Respiratory Syncytial virus type B Glycoprotein antibody [HL1733] (NBP3-25633) diluted at 1:1000. The HRP-conjugated anti-rabbit IgG antibody was used to detect the primary antibody.



Immunocytochemistry/Immunofluorescence: Respiratory Syncytial Virus Glycoprotein G Antibody (HL1733) - Azide and BSA Free [NBP3-25633] - Respiratory Syncytial virus type B Glycoprotein antibody [HL1733] detects Respiratory Syncytial virus type B Glycoprotein protein by immunofluorescent analysis. Sample: Mock and transfected 293T cells were fixed in 4% paraformaldehyde at RT for 15 min. Green: Respiratory Syncytial virus type B Glycoprotein stained by Respiratory Syncytial virus type B Glycoprotein antibody [HL1733] (NBP3-25633) diluted at 1:500. Blue: Fluoroshield with DAPI .



Immunohistochemistry-Paraffin: Respiratory Syncytial Virus Glycoprotein G Antibody (HL1733) - Azide and BSA Free [NBP3-25633] - Respiratory Syncytial virus type B G protein antibody [HL1733] detects Respiratory Syncytial virus type B G protein protein by immunohistochemical analysis. Sample:Paraffin-embedded mock and RSV Type B Major surface glycoprotein G transfected 293T. Green: Respiratory Syncytial virus type B G protein stained by Respiratory Syncytial virus type B G protein antibody [HL1733] (NBP3-25633) diluted at 1:1000. Blue: Fluoroshield with DAPI . Antigen Retrieval: Citrate buffer, pH 6.0, 15 min





Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: nb-technical@bio-
techne.com
Orders: nb-customerservice@bio-techne.com
General: novus@novusbio.com

Products Related to NBP3-25633

HAF008	Goat anti-Rabbit IgG Secondary Antibody [HRP]
NB7160	Goat anti-Rabbit IgG (H+L) Secondary Antibody [HRP]
NBP2-24891	Rabbit IgG Isotype Control

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/NBP3-25633

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

