

Product Datasheet

SMYD3 Antibody (HL2464) - Azide and BSA Free NBP3-25688

Unit Size: 100 ul

Store at 4C short term. Aliquot and store at -20C long term. Avoid freeze-thaw cycles.

www.novusbio.com



technical@novusbio.com

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/NBP3-25688

Updated 9/9/2025 v.20.1

Earn rewards for product
reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/NBP3-25688



NBP3-25688

SMYD3 Antibody (HL2464) - Azide and BSA Free


Product Information	
Unit Size	100 ul
Concentration	Concentrations vary lot to lot. See vial label for concentration. If unlisted please contact technical services.
Storage	Store at 4C short term. Aliquot and store at -20C long term. Avoid freeze-thaw cycles.
Clonality	Monoclonal
Clone	HL2464
Preservative	No Preservative
Isotype	IgG
Purity	Protein A purified
Buffer	PBS

Product Description	
Description	Novus Biologicals Rabbit SMYD3 Antibody (HL2464) - Azide and BSA Free (NBP3-25688) is a recombinant monoclonal antibody validated for use in IHC and WB. All Novus Biologicals antibodies are covered by our 100% guarantee.
Host	Rabbit
Gene ID	64754
Gene Symbol	SMYD3
Species	Human, Mouse, Rat, Zebrafish
Immunogen	Recombinant fragment of human SMYD3.


Product Application Details	
Applications	Western Blot, Immunohistochemistry-Paraffin, Immunohistochemistry
Recommended Dilutions	Western Blot 1:500-1:3000, Immunohistochemistry, Immunohistochemistry-Paraffin

Images


Western Blot: SMYD3 Antibody (HL2464) - Azide and BSA Free [NBP3-25688] - Whole cell extract (30 ug) was separated by 10% SDS-PAGE, and the membrane was blotted with SMYD3 antibody [HL2464] (NBP3-25688) diluted at 1:1000. The HRP-conjugated anti-rabbit IgG antibody was used to detect the primary antibody, and the signal was developed with Trident ECL plus-Enhanced.




Western Blot: SMYD3 Antibody (HL2464) - Azide and BSA Free [NBP3-25688] - Various whole cell extracts (30 ug) were separated by 10% SDS-PAGE, and the membrane was blotted with SMYD3 antibody [HL2464] (NBP3-25688) diluted at 1:1000. The HRP-conjugated anti-rabbit IgG antibody was used to detect the primary antibody. Corresponding RNA expression data for the same cell lines are based on Human Protein Atlas program.




Immunohistochemistry-Paraffin: SMYD3 Antibody (HL2464) - Azide and BSA Free [NBP3-25688] - SMYD3 antibody [HL2464] detects SMYD3 protein by immunohistochemical analysis. Sample: Paraffin-embedded rat tissues. SMYD3 stained by SMYD3 antibody [HL2464] (NBP3-25688) diluted at 1:100. Antigen Retrieval: Citrate buffer, pH 6.0, 15 min



Immunohistochemistry-Paraffin: SMYD3 Antibody (HL2464) - Azide and BSA Free [NBP3-25688] - SMYD3 antibody [HL2464] detects SMYD3 protein at cytoplasm by immunohistochemical analysis. Sample: Paraffin-embedded mouse colon. SMYD3 stained by SMYD3 antibody [HL2464] (NBP3-25688) diluted at 1:100. Antigen Retrieval: Citrate buffer, pH 6.0, 15 min



Western Blot: SMYD3 Antibody (HL2464) - Azide and BSA Free [NBP3-25688] - Whole zebrafish extract (50 ug) was separated by 10% SDS-PAGE, and the membrane was blotted with SMYD3 antibody [HL2464] (NBP3-25688) diluted at 1:1000. The HRP-conjugated anti-rabbit IgG antibody was used to detect the primary antibody, and the signal was developed with Trident femto Western HRP Substrate.





Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: nb-technical@bio-techne.com
Orders: nb-customerservice@bio-techne.com
General: novus@novusbio.com

Products Related to NBP3-25688

HAF008	Goat anti-Rabbit IgG Secondary Antibody [HRP]
NB7160	Goat anti-Rabbit IgG (H+L) Secondary Antibody [HRP]
NBP2-24891	Rabbit IgG Isotype Control
H00064754-P01-10ug	Recombinant Human SMYD3 GST (N-Term) Protein

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/NBP3-25688

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

