

Product Datasheet

Histone H4 [ac Lys20] Antibody (RM205) [HRP] NBP3-26013H

Unit Size: 0.1 ml

Store at 4C in the dark.

www.novusbio.com



technical@novusbio.com

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/NBP3-26013H

Updated 2/7/2024 v.20.1

Earn rewards for product
reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/NBP3-26013H



NBP3-26013H

Histone H4 [ac Lys20] Antibody (RM205) [HRP]

| Product Information | |
|------------------------------------|------------------------------------------------------------------------------------------------------------------------------|
| Unit Size | 0.1 ml |
| Concentration | Please see the vial label for concentration. If unlisted please contact technical services. |
| Storage | Store at 4C in the dark. |
| Clonality | Monoclonal |
| Clone | RM205 |
| Preservative | No Preservative |
| Isotype | IgG |
| Conjugate | HRP |
| Purity | Protein A purified |
| Buffer | PBS |
| Product Description | |
| Host | Rabbit |
| Gene ID | 121504 |
| Gene Symbol | H4C16 |
| Species | Human, Vertebrate |
| Specificity/Sensitivity | This antibody reacts to Histone H4 acetylated at Lysine 20. No cross reactivity with other acetylated Lysines in Histone H4. |
| Immunogen | An acetyl-peptide corresponding to the Histone H4 [ac Lys20] |
| Product Application Details | |
| Applications | Western Blot, ELISA, Immunocytochemistry/ Immunofluorescence, Chromatin Immunoprecipitation (ChIP), Multiplex Immunoassay |
| Recommended Dilutions | Western Blot, ELISA, Immunocytochemistry/ Immunofluorescence, Chromatin Immunoprecipitation (ChIP), Multiplex Immunoassay |
| Application Notes | Optimal dilution of this antibody should be experimentally determined. |



Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: nb-technical@bio-techne.com
Orders: nb-customerservice@bio-techne.com
General: novus@novusbio.com

Products Related to NBP3-26013H

| | |
|--------------|--------------------------------------------------------|
| NBP2-24891H | Rabbit IgG Isotype Control [HRP] |
| NB21-2011PEP | Histone H4 [Monomethyl Arg3] Antibody Blocking Peptide |
| 2695-SE-010 | Airway Trypsin-like Protease/HAT/TMPRSS11D |
| NB21-2089PEP | Histone H4 [Dimethyl Lys20] Antibody Blocking Peptide |

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/NBP3-26013H

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

