

# Product Datasheet

## IL-12/IL-23 p40 Antibody (171) [HRP] NBP3-26107H

Unit Size: 0.1 ml

Store at 4C in the dark.

[www.novusbio.com](http://www.novusbio.com)



[technical@novusbio.com](mailto:technical@novusbio.com)

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:  
[www.novusbio.com/NBP3-26107H](http://www.novusbio.com/NBP3-26107H)

Updated 2/15/2024 v.20.1

Earn rewards for product  
reviews and publications.

Submit a publication at [www.novusbio.com/publications](http://www.novusbio.com/publications)

Submit a review at [www.novusbio.com/reviews/destination/NBP3-26107H](http://www.novusbio.com/reviews/destination/NBP3-26107H)



**NBP3-26107H**

IL-12/IL-23 p40 Antibody (171) [HRP]

| <b>Product Information</b> |                                                                                             |
|----------------------------|---------------------------------------------------------------------------------------------|
| <b>Unit Size</b>           | 0.1 ml                                                                                      |
| <b>Concentration</b>       | Please see the vial label for concentration. If unlisted please contact technical services. |
| <b>Storage</b>             | Store at 4C in the dark.                                                                    |
| <b>Clonality</b>           | Monoclonal                                                                                  |
| <b>Clone</b>               | 171                                                                                         |
| <b>Preservative</b>        | No Preservative                                                                             |
| <b>Isotype</b>             | IgG                                                                                         |
| <b>Conjugate</b>           | HRP                                                                                         |
| <b>Purity</b>              | Protein A purified                                                                          |
| <b>Buffer</b>              | PBS                                                                                         |

| <b>Product Description</b> |                                                                                                                                      |
|----------------------------|--------------------------------------------------------------------------------------------------------------------------------------|
| <b>Host</b>                | Rabbit                                                                                                                               |
| <b>Gene ID</b>             | 3593                                                                                                                                 |
| <b>Gene Symbol</b>         | IL12B                                                                                                                                |
| <b>Species</b>             | Mouse                                                                                                                                |
| <b>Immunogen</b>           | This antibody was obtained from a rabbit immunized with purified, recombinant Mouse IL-12 p40/IL12B (Uniprot#: P43432; Met1-Ser335). |

| <b>Product Application Details</b> |                                                                        |
|------------------------------------|------------------------------------------------------------------------|
| <b>Applications</b>                | ELISA                                                                  |
| <b>Recommended Dilutions</b>       | ELISA                                                                  |
| <b>Application Notes</b>           | Optimal dilution of this antibody should be experimentally determined. |





### **Novus Biologicals USA**

10730 E. Briarwood Avenue  
Centennial, CO 80112  
USA  
Phone: 303.730.1950  
Toll Free: 1.888.506.6887  
Fax: 303.730.1966  
nb-customerservice@bio-techne.com

### **Bio-Techne Canada**

21 Canmotor Ave  
Toronto, ON M8Z 4E6  
Canada  
Phone: 905.827.6400  
Toll Free: 855.668.8722  
Fax: 905.827.6402  
canada.inquires@bio-techne.com

### **Bio-Techne Ltd**

19 Barton Lane  
Abingdon Science Park  
Abingdon, OX14 3NB, United Kingdom  
Phone: (44) (0) 1235 529449  
Free Phone: 0800 37 34 15  
Fax: (44) (0) 1235 533420  
info.EMEA@bio-techne.com

### **General Contact Information**

www.novusbio.com  
Technical Support: nb-technical@bio-techne.com  
Orders: nb-customerservice@bio-techne.com  
General: novus@novusbio.com

### **Products Related to NBP3-26107H**

---

|             |                                  |
|-------------|----------------------------------|
| NBP2-24891H | Rabbit IgG Isotype Control [HRP] |
| 210-TA-005  | TNF-alpha [Unconjugated]         |
| DY499-05    | IL-12/IL-23 p40 [Biotin]         |
| M6000B-1    | IL-6 [HRP]                       |

---

### **Limitations**

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit [www.novusbio.com/guarantee](http://www.novusbio.com/guarantee)

Earn gift cards/discounts by submitting a review: [www.novusbio.com/reviews/submit/NBP3-26107H](http://www.novusbio.com/reviews/submit/NBP3-26107H)

Earn gift cards/discounts by submitting a publication using this product:  
[www.novusbio.com/publications](http://www.novusbio.com/publications)

