

# Product Datasheet

## PAUF/ZG16B Antibody (6E5) - BSA Free NBP3-26244-100ul

Unit Size: 100 ul

Store at -20 to -70C. Avoid freeze-thaw cycles.

[www.novusbio.com](http://www.novusbio.com)



[technical@novusbio.com](mailto:technical@novusbio.com)

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:  
[www.novusbio.com/NBP3-26244](http://www.novusbio.com/NBP3-26244)

Updated 2/24/2026 v.20.1

Earn rewards for product  
reviews and publications.

Submit a publication at [www.novusbio.com/publications](http://www.novusbio.com/publications)

Submit a review at [www.novusbio.com/reviews/destination/NBP3-26244](http://www.novusbio.com/reviews/destination/NBP3-26244)



**NBP3-26244-100ul**

PAUF/ZG16B Antibody (6E5) - BSA Free

Product Information	
Unit Size	100 ul
Concentration	Please see the vial label for concentration. If unlisted please contact technical services.
Storage	Store at -20 to -70C. Avoid freeze-thaw cycles.
Clonality	Monoclonal
Clone	6E5
Preservative	0.03% Proclin 300
Isotype	IgG1
Purity	Affinity purified
Buffer	50% Glycerol, 0.01M PBS, PH 7.4


Product Description	
Description	Novus Biologicals Human PAUF/ZG16B Antibody (6E5) - BSA Free (NBP3-26244) is a recombinant monoclonal antibody validated for use in ELISA. All Novus Biologicals antibodies are covered by our 100% guarantee.
Host	Human
Gene ID	124220
Gene Symbol	ZG16B
Species	Human, Cynomolgus Monkey
Immunogen	Recombinant Human PAUF/ZG16B protein [UniProt Q96DA0]

Product Application Details	
Applications	ELISA
Recommended Dilutions	ELISA




## Images

ELISA: PAUF/ZG16B Antibody (6E5) [NBP3-26244] - The Binding Activity of Human PAUF/ZG16B with PAUF/ZG16B Antibody.  
 Activity: Measured by its binding ability in a functional ELISA.  
 Immobilized Human PAUF/ZG16B at 2  $\mu\text{g}/\text{mL}$  can bind PAUF/ZG16B Antibody, the EC50 is 24.13-46.04  $\text{ng}/\text{mL}$ .



ELISA: PAUF/ZG16B Antibody (6E5) [NBP3-26244] -The Binding Activity of Macaca fascicularis ZG16B with Anti-ZG16B recombinant antibody  
 Activity: Measured by its binding ability in a functional ELISA.  
 Immobilized Macaca fascicularis ZG16B at 2  $\mu\text{g}/\text{mL}$  can bind Anti-ZG16B recombinant antibody , the EC50 is 20.15-28.98  $\text{ng}/\text{mL}$ .





### **Novus Biologicals USA**

10730 E. Briarwood Avenue  
Centennial, CO 80112  
USA  
Phone: 303.730.1950  
Toll Free: 1.888.506.6887  
Fax: 303.730.1966  
nb-customerservice@bio-techne.com

### **Bio-Techne Canada**

21 Canmotor Ave  
Toronto, ON M8Z 4E6  
Canada  
Phone: 905.827.6400  
Toll Free: 855.668.8722  
Fax: 905.827.6402  
canada.inquires@bio-techne.com

### **Bio-Techne Ltd**

19 Barton Lane  
Abingdon Science Park  
Abingdon, OX14 3NB, United Kingdom  
Phone: (44) (0) 1235 529449  
Free Phone: 0800 37 34 15  
Fax: (44) (0) 1235 533420  
info.EMEA@bio-techne.com

### **General Contact Information**

www.novusbio.com  
Technical Support: nb-technical@bio-techne.com  
Orders: nb-customerservice@bio-techne.com  
General: novus@novusbio.com

### **Products Related to NBP3-26244-100ul**

---

G-102-C	Goat anti-Human IgG Secondary Antibody [Unconjugated]
NB7446	Goat anti-Human IgG Fc Secondary Antibody
DDXCH01P-100	Human IgG1 Isotype Control
H00124220-P01-10ug	Recombinant Human PAUF/ZG16B GST (N-Term) Protein

---

### **Limitations**

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit [www.novusbio.com/guarantee](http://www.novusbio.com/guarantee)

Earn gift cards/discounts by submitting a review: [www.novusbio.com/reviews/submit/NBP3-26244](http://www.novusbio.com/reviews/submit/NBP3-26244)

Earn gift cards/discounts by submitting a publication using this product:  
[www.novusbio.com/publications](http://www.novusbio.com/publications)

