

Product Datasheet

Fc gamma RIIIA/CD16a Antibody (28C11) NBP3-26586-100ul

Unit Size: 100 ul

Store at -20 to -70C. Avoid freeze-thaw cycles.

www.novusbio.com



technical@novusbio.com

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/NBP3-26586

Updated 2/21/2024 v.20.1

Earn rewards for product
reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/NBP3-26586



NBP3-26586-100ul

Fc gamma RIIIA/CD16a Antibody (28C11)


| Product Information | |
|---------------------|---|
| Unit Size | 100 ul |
| Concentration | Please see the vial label for concentration. If unlisted please contact technical services. |
| Storage | Store at -20 to -70C. Avoid freeze-thaw cycles. |
| Clonality | Monoclonal |
| Clone | 28C11 |
| Preservative | 0.03% Proclin 300 |
| Isotype | IgG |
| Purity | Affinity purified |
| Buffer | 50% Glycerol, 0.01M PBS, PH 7.4 |

| Product Description | |
|---------------------|---|
| Host | Rabbit |
| Gene ID | 2214 |
| Gene Symbol | FCGR3A |
| Species | Human |
| Immunogen | Recombinant Human Fc gamma RIIIA/CD16a protein [UniProt P08637] |

| Product Application Details | |
|-----------------------------|---|
| Applications | ELISA, Immunocytochemistry/ Immunofluorescence |
| Recommended Dilutions | ELISA, Immunocytochemistry/ Immunofluorescence 1:50-1:200 |

Images

Immunocytochemistry/Immunofluorescence: Fc gamma RIIIA/CD16a Antibody (28C11) [NBP3-26586] - Staining of Hela cell with Fc gamma RIIIA/CD16a Antibody (28C11) at 1:30, counter-stained with DAPI. The cells were fixed in 4% formaldehyde and blocked in 10% normal Goat Serum. The cells were then incubated with the antibody overnight at 4C. The secondary antibody was Alexa Fluor 488-conjugated Goat Anti-Rabbit IgG (H+L).





Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: nb-technical@bio-techne.com
Orders: nb-customerservice@bio-techne.com
General: novus@novusbio.com

Products Related to NBP3-26586-100ul

| | |
|---------------|---|
| HAF008 | Goat anti-Rabbit IgG Secondary Antibody [HRP] |
| NB7160 | Goat anti-Rabbit IgG (H+L) Secondary Antibody [HRP] |
| NBP2-24891 | Rabbit IgG Isotype Control |
| NBP2-33753PEP | Fc gamma RIIIA/CD16a Recombinant Protein Antigen |

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/NBP3-26586

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

