

Product Datasheet

Cdk6 Antibody (8G3)

NBP3-26663-100ul

Unit Size: 100 ul

Store at -20 to -70C. Avoid freeze-thaw cycles.

www.novusbio.com



technical@novusbio.com

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/NBP3-26663

Updated 1/7/2025 v.20.1

Earn rewards for product
reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/NBP3-26663




NBP3-26663-100ul

Cdk6 Antibody (8G3)


Product Information	
Unit Size	100 ul
Concentration	Please see the vial label for concentration. If unlisted please contact technical services.
Storage	Store at -20 to -70C. Avoid freeze-thaw cycles.
Clonality	Monoclonal
Clone	8G3
Preservative	0.02% Sodium Azide
Isotype	IgG
Purity	Affinity purified
Buffer	PBS, pH 7.4, 150mM NaCl, and 50% glycerol
Product Description	
Host	Rabbit
Gene ID	1021
Gene Symbol	CDK6
Species	Human
Immunogen	A synthesized peptide derived from Human Cdk6 [UniProt Q00534]
Product Application Details	
Applications	Western Blot, ELISA, Immunoprecipitation
Recommended Dilutions	Western Blot 1:500-1:5000, ELISA, Immunoprecipitation 1:200-1:1000

Images

Immunoprecipitation: Cdk6 Antibody (8G3) [NBP3-26663] - Immunoprecipitating Cdk6 in K562 whole cell lysate.
 Lane 1: Rabbit control IgG instead of NBP3-26663 in K562 whole cell lysate. For western blotting, a HRP-conjugated Protein G antibody was used as the secondary antibody (1/2000).
 Lane 2: NBP3-26663 (2ug)+ K562 whole cell lysate (500ug).
 Lane 3: K562 whole cell lysate (10ug).



Western Blot: Cdk6 Antibody (8G3) [NBP3-26663] - Positive Western Blot detected in: Hela whole cell lysate, Jurkat whole cell lysate, K562 whole cell lysate, 293 whole cell lysate.
 All lanes: Cdk6 Antibody at 1: 2000
 Secondary: Goat polyclonal to rabbit IgG at 1/50000 dilution.
 Predicted band size: 37 kDa
 Observed band size: 37 kDa





Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: nb-technical@bio-techne.com
Orders: nb-customerservice@bio-techne.com
General: novus@novusbio.com

Products Related to NBP3-26663-100ul

HAF008	Goat anti-Rabbit IgG Secondary Antibody [HRP]
NB7160	Goat anti-Rabbit IgG (H+L) Secondary Antibody [HRP]
NBP2-24891	Rabbit IgG Isotype Control
H00001021-P01-10ug	Recombinant Human Cdk6 GST (N-Term) Protein

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/NBP3-26663

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

