

# Product Datasheet

## UIP1 Antibody (SR4274)

### NBP3-27830-100ul

Unit Size: 100 ul

Store at 4C short term. Aliquot and store at -20C long term. Avoid freeze-thaw cycles.

[www.novusbio.com](http://www.novusbio.com)



[technical@novusbio.com](mailto:technical@novusbio.com)

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:  
[www.novusbio.com/NBP3-27830](http://www.novusbio.com/NBP3-27830)

Updated 5/5/2024 v.20.1

Earn rewards for product  
reviews and publications.

Submit a publication at [www.novusbio.com/publications](http://www.novusbio.com/publications)

Submit a review at [www.novusbio.com/reviews/destination/NBP3-27830](http://www.novusbio.com/reviews/destination/NBP3-27830)



**NBP3-27830-100ul**

UIP1 Antibody (SR4274)

| Product Information         |                                                                                             |
|-----------------------------|---------------------------------------------------------------------------------------------|
| Unit Size                   | 100 ul                                                                                      |
| Concentration               | Please see the vial label for concentration. If unlisted please contact technical services. |
| Storage                     | Store at 4C short term. Aliquot and store at -20C long term. Avoid freeze-thaw cycles.      |
| Clonality                   | Monoclonal                                                                                  |
| Clone                       | SR4274                                                                                      |
| Preservative                | 0.02% Sodium Azide                                                                          |
| Isotype                     | IgG                                                                                         |
| Purity                      | Affinity purified                                                                           |
| Buffer                      | PBS (pH 7.4), 150mM NaCl, 50% glycerol                                                      |
| Target Molecular Weight     | 37 kDa                                                                                      |
| Product Description         |                                                                                             |
| Host                        | Rabbit                                                                                      |
| Gene Symbol                 | HAUS7                                                                                       |
| Species                     | Human                                                                                       |
| Specificity/Sensitivity     | Detects endogenous levels of total UIP1                                                     |
| Immunogen                   | A synthesized peptide derived from human UIP1 (Uniprot #: Q99871)                           |
| Product Application Details |                                                                                             |
| Applications                | Western Blot, Immunohistochemistry, Immunoprecipitation                                     |
| Recommended Dilutions       | Western Blot 1:500-1:2000, Immunohistochemistry 1:50-1:200, Immunoprecipitation 1:50        |





### **Novus Biologicals USA**

10730 E. Briarwood Avenue  
Centennial, CO 80112  
USA  
Phone: 303.730.1950  
Toll Free: 1.888.506.6887  
Fax: 303.730.1966  
nb-customerservice@bio-techne.com

### **Bio-Techne Canada**

21 Canmotor Ave  
Toronto, ON M8Z 4E6  
Canada  
Phone: 905.827.6400  
Toll Free: 855.668.8722  
Fax: 905.827.6402  
canada.inquires@bio-techne.com

### **Bio-Techne Ltd**

19 Barton Lane  
Abingdon Science Park  
Abingdon, OX14 3NB, United Kingdom  
Phone: (44) (0) 1235 529449  
Free Phone: 0800 37 34 15  
Fax: (44) (0) 1235 533420  
info.EMEA@bio-techne.com

### **General Contact Information**

[www.novusbio.com](http://www.novusbio.com)  
Technical Support: nb-technical@bio-techne.com  
Orders: nb-customerservice@bio-techne.com  
General: novus@novusbio.com

### **Products Related to NBP3-27830-100ul**

---

|                    |                                                     |
|--------------------|-----------------------------------------------------|
| HAF008             | Goat anti-Rabbit IgG Secondary Antibody [HRP]       |
| NB7160             | Goat anti-Rabbit IgG (H+L) Secondary Antibody [HRP] |
| NBP2-24891         | Rabbit IgG Isotype Control                          |
| H00055559-P01-10ug | Recombinant Human UIP1 GST (N-Term) Protein         |

---

### **Limitations**

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit [www.novusbio.com/guarantee](http://www.novusbio.com/guarantee)

Earn gift cards/discounts by submitting a review: [www.novusbio.com/reviews/submit/NBP3-27830](http://www.novusbio.com/reviews/submit/NBP3-27830)

Earn gift cards/discounts by submitting a publication using this product:  
[www.novusbio.com/publications](http://www.novusbio.com/publications)

