

Product Datasheet

IFIT3 Antibody - BSA Free NBP3-27893

Unit Size: 0.1 mg

Store at 4C short term. Aliquot and store at -20C long term. Avoid freeze-thaw cycles.

www.novusbio.com



technical@novusbio.com

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/NBP3-27893

Updated 9/9/2025 v.20.1

Earn rewards for product
reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/NBP3-27893



NBP3-27893

IFIT3 Antibody - BSA Free

Product Information

Unit Size	0.1 mg
Concentration	1 mg/ml
Storage	Store at 4C short term. Aliquot and store at -20C long term. Avoid freeze-thaw cycles.
Clonality	Polyclonal
Preservative	0.02% Sodium Azide
Purity	Peptide affinity purified
Buffer	PBS

Product Description

Description	Novus Biologicals Rabbit IFIT3 Antibody - BSA Free (NBP3-27893) is a polyclonal antibody validated for use in IHC, WB and ELISA. All Novus Biologicals antibodies are covered by our 100% guarantee.
Host	Rabbit
Gene ID	3437
Gene Symbol	IFIT3
Species	Human
Specificity/Sensitivity	IFIT3 antibody is human specific. This antibody is predicted to not cross-react with other members of the IFIT protein family.
Immunogen	IFIT3 antibody was raised against an 18 amino acid peptide near the carboxy terminus of human IFIT3. The immunogen is located within amino acids 380 - 430 of IFIT3 .


Product Application Details

Applications	Western Blot, Immunohistochemistry-Paraffin, ELISA
Recommended Dilutions	Western Blot 1 - 2 ug/mL, ELISA, Immunohistochemistry-Paraffin 5 ug/mL




Images

Western Blot: IFIT3 Antibody [NBP3-27893] - Western blot analysis of IFIT3 in HepG2 cell lysate with IFIT3 antibody at 1 u/ml.



Immunohistochemistry: IFIT3 Antibody [NBP3-27893] - Immunohistochemistry of IFIT3 in human liver tissue with IFIT3 antibody at 5 u/mL.





Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: nb-technical@bio-techne.com
Orders: nb-customerservice@bio-techne.com
General: novus@novusbio.com

Products Related to NBP3-27893

HAF008	Goat anti-Rabbit IgG Secondary Antibody [HRP]
NB7160	Goat anti-Rabbit IgG (H+L) Secondary Antibody [HRP]
NBP2-52117-0.1mg	Recombinant Human IFIT3 His Protein
210-TA-005	TNF-alpha [Unconjugated]

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/NBP3-27893

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

