

# Product Datasheet

## IFIT3 Antibody [Janelia Fluor® 646] NBP3-27893JF646

Unit Size: 0.1 ml

Store at 4C in the dark.

[www.novusbio.com](http://www.novusbio.com)



[technical@novusbio.com](mailto:technical@novusbio.com)

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:  
[www.novusbio.com/NBP3-27893JF646](http://www.novusbio.com/NBP3-27893JF646)

Updated 3/24/2025 v.20.1

Earn rewards for product  
reviews and publications.

Submit a publication at [www.novusbio.com/publications](http://www.novusbio.com/publications)

Submit a review at [www.novusbio.com/reviews/destination/NBP3-27893JF646](http://www.novusbio.com/reviews/destination/NBP3-27893JF646)



**NBP3-27893JF646**

IFIT3 Antibody [Janelia Fluor® 646]

| Product Information |                                                                                             |
|---------------------|---------------------------------------------------------------------------------------------|
| Unit Size           | 0.1 ml                                                                                      |
| Concentration       | Please see the vial label for concentration. If unlisted please contact technical services. |
| Storage             | Store at 4C in the dark.                                                                    |
| Clonality           | Polyclonal                                                                                  |
| Preservative        | 0.05% Sodium Azide                                                                          |
| Conjugate           | Janelia Fluor 646                                                                           |
| Purity              | Peptide affinity purified                                                                   |
| Buffer              | 50mM Sodium Borate                                                                          |

| Product Description     |                                                                                                                                                                       |
|-------------------------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Host                    | Rabbit                                                                                                                                                                |
| Gene ID                 | 3437                                                                                                                                                                  |
| Gene Symbol             | IFIT3                                                                                                                                                                 |
| Species                 | Human                                                                                                                                                                 |
| Specificity/Sensitivity | IFIT3 antibody is human specific. This antibody is predicted to not cross-react with other members of the IFIT protein family.                                        |
| Immunogen               | IFIT3 antibody was raised against an 18 amino acid peptide near the carboxy terminus of human IFIT3. The immunogen is located within amino acids 380 - 430 of IFIT3 . |
| Notes                   | Sold under license from the Howard Hughes Medical Institute, Janelia Research Campus.                                                                                 |

| Product Application Details |                                                                        |
|-----------------------------|------------------------------------------------------------------------|
| Applications                | Western Blot, ELISA, Immunohistochemistry-Paraffin                     |
| Recommended Dilutions       | Western Blot, ELISA, Immunohistochemistry-Paraffin                     |
| Application Notes           | Optimal dilution of this antibody should be experimentally determined. |





### **Novus Biologicals USA**

10730 E. Briarwood Avenue  
Centennial, CO 80112  
USA  
Phone: 303.730.1950  
Toll Free: 1.888.506.6887  
Fax: 303.730.1966  
nb-customerservice@bio-techne.com

### **Bio-Techne Canada**

21 Canmotor Ave  
Toronto, ON M8Z 4E6  
Canada  
Phone: 905.827.6400  
Toll Free: 855.668.8722  
Fax: 905.827.6402  
canada.inquires@bio-techne.com

### **Bio-Techne Ltd**

19 Barton Lane  
Abingdon Science Park  
Abingdon, OX14 3NB, United Kingdom  
Phone: (44) (0) 1235 529449  
Free Phone: 0800 37 34 15  
Fax: (44) (0) 1235 533420  
info.EMEA@bio-techne.com

### **General Contact Information**

[www.novusbio.com](http://www.novusbio.com)  
Technical Support: nb-technical@bio-techne.com  
Orders: nb-customerservice@bio-techne.com  
General: novus@novusbio.com

### **Products Related to NBP3-27893JF646**

---

|                  |                                     |
|------------------|-------------------------------------|
| NBP2-52117-0.1mg | Recombinant Human IFIT3 His Protein |
| 210-TA-005       | TNF-alpha [Unconjugated]            |
| NBL1-11837       | IFIT3 Overexpression Lysate         |
| NB200-106        | p19ARF/CDKN2A Antibody - BSA Free   |

---

### **Limitations**

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit [www.novusbio.com/guarantee](http://www.novusbio.com/guarantee)

Earn gift cards/discounts by submitting a review: [www.novusbio.com/reviews/submit/NBP3-27893JF646](http://www.novusbio.com/reviews/submit/NBP3-27893JF646)

Earn gift cards/discounts by submitting a publication using this product:  
[www.novusbio.com/publications](http://www.novusbio.com/publications)

