

# Product Datasheet

## **B7-H6 Antibody - Humanized - Low Endotoxin, Azide and BSA Free**

### **NBP3-27971-100ug**

Unit Size: 100 ug

Store at -20C in powder form. Store at -80C once reconstituted.

[www.novusbio.com](http://www.novusbio.com)



[technical@novusbio.com](mailto:technical@novusbio.com)

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:  
[www.novusbio.com/NBP3-27971](http://www.novusbio.com/NBP3-27971)

Updated 12/8/2025 v.20.1

Earn rewards for product reviews and publications.

Submit a publication at [www.novusbio.com/publications](http://www.novusbio.com/publications)

Submit a review at [www.novusbio.com/reviews/destination/NBP3-27971](http://www.novusbio.com/reviews/destination/NBP3-27971)



**NBP3-27971-100ug**

B7-H6 Antibody - Humanized - Low Endotoxin, Azide and BSA Free

**Product Information**

|                                    |                                                                                                                                       |
|------------------------------------|---------------------------------------------------------------------------------------------------------------------------------------|
| <b>Unit Size</b>                   | 100 ug                                                                                                                                |
| <b>Concentration</b>               | LYOPH mg/ml                                                                                                                           |
| <b>Storage</b>                     | Store at -20C in powder form. Store at -80C once reconstituted.                                                                       |
| <b>Clonality</b>                   | Monoclonal                                                                                                                            |
| <b>Preservative</b>                | No Preservative                                                                                                                       |
| <b>Reconstitution Instructions</b> | Reconstitute with sterile, distilled water to a final concentration of 1 mg/ml. Gently shake to solubilize completely. Do not vortex. |
| <b>Isotype</b>                     | IgG1                                                                                                                                  |
| <b>Purity</b>                      | Protein A purified                                                                                                                    |
| <b>Buffer</b>                      | Lyophilized from 25mM histidine, 8% sucrose, 0.01% Tween80 (pH6.2)                                                                    |

**Product Description**

|                       |                                                                                                                                                                                                                                                                                                   |
|-----------------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Description</b>    | <p>It has a predicted MW of 145.5 kDa.</p> <p>This product is shipped at ambient temperature. Upon receipt, store immediately at -20C or lower for 12 months in a lyophilized state. - 80C for 3 months after reconstitution. Avoid repeated freeze-thaw cycles.</p> <p>Patent #: US9790278B2</p> |
| <b>Host</b>           | Human                                                                                                                                                                                                                                                                                             |
| <b>Gene ID</b>        | 374383                                                                                                                                                                                                                                                                                            |
| <b>Gene Symbol</b>    | NCR3LG1                                                                                                                                                                                                                                                                                           |
| <b>Species</b>        | Human                                                                                                                                                                                                                                                                                             |
| <b>Immunogen</b>      | B7-H6 / NCR3LG1                                                                                                                                                                                                                                                                                   |
| <b>Endotoxin Note</b> | < 0.001EU/ug,determined by LAL method.                                                                                                                                                                                                                                                            |

**Product Application Details**

|                              |                                   |
|------------------------------|-----------------------------------|
| <b>Applications</b>          | ELISA, Flow Cytometry, Functional |
| <b>Recommended Dilutions</b> | Flow Cytometry, ELISA, Functional |





### **Novus Biologicals USA**

10730 E. Briarwood Avenue  
Centennial, CO 80112  
USA  
Phone: 303.730.1950  
Toll Free: 1.888.506.6887  
Fax: 303.730.1966  
nb-customerservice@bio-techne.com

### **Bio-Techne Canada**

21 Canmotor Ave  
Toronto, ON M8Z 4E6  
Canada  
Phone: 905.827.6400  
Toll Free: 855.668.8722  
Fax: 905.827.6402  
canada.inquires@bio-techne.com

### **Bio-Techne Ltd**

19 Barton Lane  
Abingdon Science Park  
Abingdon, OX14 3NB, United Kingdom  
Phone: (44) (0) 1235 529449  
Free Phone: 0800 37 34 15  
Fax: (44) (0) 1235 533420  
info.EMEA@bio-techne.com

### **General Contact Information**

[www.novusbio.com](http://www.novusbio.com)  
Technical Support: [nb-technical@bio-techne.com](mailto:nb-technical@bio-techne.com)  
Orders: [nb-customerservice@bio-techne.com](mailto:nb-customerservice@bio-techne.com)  
General: [novus@novusbio.com](mailto:novus@novusbio.com)

### **Products Related to NBP3-27971-100ug**

---

|              |                                                       |
|--------------|-------------------------------------------------------|
| G-102-C      | Goat anti-Human IgG Secondary Antibody [Unconjugated] |
| NB7446       | Goat anti-Human IgG Fc Secondary Antibody             |
| DDXCH01P-100 | Human IgG1 Isotype Control                            |
| 7268-CT-100  | CTLA-4 [Unconjugated]                                 |

---

### **Limitations**

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit [www.novusbio.com/guarantee](http://www.novusbio.com/guarantee)

Earn gift cards/discounts by submitting a review: [www.novusbio.com/reviews/submit/NBP3-27971](http://www.novusbio.com/reviews/submit/NBP3-27971)

Earn gift cards/discounts by submitting a publication using this product:  
[www.novusbio.com/publications](http://www.novusbio.com/publications)

