

Product Datasheet

Haptoglobin Antibody (KHK patent anti-Haptoglobin) - Humanized NBP3-28166-100ug

Unit Size: 100 ug

Store at -20C in powder form. Store at -80C once reconstituted.

www.novusbio.com



technical@novusbio.com

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/NBP3-28166

Updated 7/7/2024 v.20.1

Earn rewards for product
reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/NBP3-28166



NBP3-28166-100ug

Haptoglobin Antibody (KHK patent anti-Haptoglobin) - Humanized

| Product Information | |
|-----------------------------|---|
| Unit Size | 100 ug |
| Concentration | LYOPH mg/ml |
| Storage | Store at -20C in powder form. Store at -80C once reconstituted. |
| Clonality | Monoclonal |
| Clone | KHK patent anti-Haptoglobin |
| Preservative | No Preservative |
| Reconstitution Instructions | Reconstitute with sterile, distilled water to a final concentration of 1 mg/ml. Gently shake to solubilize completely. Do not vortex. |
| Isotype | IgG1 |
| Purity | Protein A purified |
| Buffer | Lyophilized from 25mM histidine, 8% sucrose, 0.01% Tween80 (pH6.2) |

| Product Description | |
|---------------------|--|
| Description | Expressed from CHO. The heavy chain type is hulgG1, and the light chain type is hukappa. It has a predicted MW of 145.5 kDa. Upon receipt, store immediately at -20C or lower for 24 months in a lyophilized state. - 80C for 3 months after reconstitution. Avoid repeated freeze-thaw cycles. |
| Host | Human |
| Gene ID | 3240 |
| Gene Symbol | HP |
| Species | Human |
| Immunogen | Haptoglobin |
| Endotoxin Note | < 0.001EU/ug,determined by LAL method. |

| Product Application Details | |
|-----------------------------|-----------------------------------|
| Applications | ELISA, Flow Cytometry, Functional |
| Recommended Dilutions | Flow Cytometry, ELISA, Functional |





Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: nb-technical@bio-techne.com
Orders: nb-customerservice@bio-techne.com
General: novus@novusbio.com

Products Related to NBP3-28166-100ug

| | |
|-------------------|---|
| G-102-C | Goat anti-Human IgG Secondary Antibody [Unconjugated] |
| NB7446 | Goat anti-Human IgG Fc Secondary Antibody |
| DDXCH01P-100 | Human IgG1 Isotype Control |
| NBP2-51611-0.05mg | Recombinant Human Haptoglobin His Protein |

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/NBP3-28166

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

