

# Product Datasheet

**alpha-Synuclein Antibody - Humanized, (Research Grade  
cinpanemab Biosimilar) - Low Endotoxin, Azide and BSA  
Free**

**NBP3-28298-100ug**

Unit Size: 100 ug

Store at 4C.

[www.novusbio.com](http://www.novusbio.com)



[technical@novusbio.com](mailto:technical@novusbio.com)

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:


[www.novusbio.com/NBP3-28298](http://www.novusbio.com/NBP3-28298)

Updated 3/13/2025 v.20.1

Earn rewards for product  
reviews and publications.

Submit a publication at [www.novusbio.com/publications](http://www.novusbio.com/publications)

Submit a review at [www.novusbio.com/reviews/destination/NBP3-28298](http://www.novusbio.com/reviews/destination/NBP3-28298)



**NBP3-28298-100ug**

alpha-Synuclein Antibody - Humanized, (Research Grade cinpanemab Biosimilar) - Low Endotoxin, Azide and BSA Free

**Product Information**

<b>Unit Size</b>	100 ug
<b>Concentration</b>	LYOPH mg/ml
<b>Storage</b>	Store at 4C.
<b>Clonality</b>	Monoclonal
<b>Preservative</b>	No Preservative
<b>Reconstitution Instructions</b>	Reconstitute with sterile, distilled water to a final concentration of 1 mg/ml. Gently shake to solubilize completely. Do not vortex.
<b>Isotype</b>	IgG1
<b>Purity</b>	Protein A purified
<b>Buffer</b>	Lyophilized from 25mM histidine, 8% sucrose, 0.01% Tween80 (pH6.2)

**Product Description**


<b>Description</b>	Expressed from CHO. The heavy chain type is hulG1, and the light chain type is hulambda. It has a predicted MW of 146.98 kDa.  Also known as 'cinpanemab'.  Upon receipt, store immediately at -20C or lower for 24 months in a lyophilized state. - 80C for 3 months after reconstitution. Avoid repeated freeze-thaw cycles.
<b>Host</b>	Human
<b>Gene ID</b>	6622
<b>Gene Symbol</b>	SNCA
<b>Species</b>	Human
<b>Immunogen</b>	Alpha-synuclein
<b>Endotoxin Note</b>	< 0.001EU/ug,determined by LAL method.

**Product Application Details**

<b>Applications</b>	ELISA, Flow Cytometry, Functional
<b>Recommended Dilutions</b>	Flow Cytometry, ELISA, Functional

**Images**

Immobilized Human Alpha-synuclein, Fc Tag at 2 µg/mL can bind alpha-Synuclein Antibody (cinpanemab) - Humanized, The EC50 was approximately 0.02682 ug/mL.





### **Novus Biologicals USA**

10730 E. Briarwood Avenue  
Centennial, CO 80112  
USA  
Phone: 303.730.1950  
Toll Free: 1.888.506.6887  
Fax: 303.730.1966  
nb-customerservice@bio-techne.com

### **Bio-Techne Canada**

21 Canmotor Ave  
Toronto, ON M8Z 4E6  
Canada  
Phone: 905.827.6400  
Toll Free: 855.668.8722  
Fax: 905.827.6402  
canada.inquires@bio-techne.com

### **Bio-Techne Ltd**

19 Barton Lane  
Abingdon Science Park  
Abingdon, OX14 3NB, United Kingdom  
Phone: (44) (0) 1235 529449  
Free Phone: 0800 37 34 15  
Fax: (44) (0) 1235 533420  
info.EMEA@bio-techne.com

### **General Contact Information**

www.novusbio.com  
Technical Support: nb-technical@bio-techne.com  
Orders: nb-customerservice@bio-techne.com  
General: novus@novusbio.com

### **Products Related to NBP3-28298-100ug**

---

G-102-C	Goat anti-Human IgG Secondary Antibody [Unconjugated]
NB7446	Goat anti-Human IgG Fc Secondary Antibody
DDXCH01P-100	Human IgG1 Isotype Control
NBP2-54789-100ug	Recombinant Human alpha-Synuclein Active, Pre-formed Fibrils, (Type 1) Protein

---

### **Limitations**

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit [www.novusbio.com/guarantee](http://www.novusbio.com/guarantee)

Earn gift cards/discounts by submitting a review: [www.novusbio.com/reviews/submit/NBP3-28298](http://www.novusbio.com/reviews/submit/NBP3-28298)

Earn gift cards/discounts by submitting a publication using this product:  
[www.novusbio.com/publications](http://www.novusbio.com/publications)

