

# Product Datasheet

## **ST2/IL-33R Antibody - IgG2SA, (Research Grade melrilimab Biosimilar) - Low Endotoxin, Azide and BSA Free NBP3-28740-100ug**

Unit Size: 100 ug

Store at 4C.

[www.novusbio.com](http://www.novusbio.com)



[technical@novusbio.com](mailto:technical@novusbio.com)


Protocols, Publications, Related Products, Reviews, Research Tools and Images at:  
[www.novusbio.com/NBP3-28740](http://www.novusbio.com/NBP3-28740)

Updated 3/13/2025 v.20.1

Earn rewards for product  
reviews and publications.

Submit a publication at [www.novusbio.com/publications](http://www.novusbio.com/publications)

Submit a review at [www.novusbio.com/reviews/destination/NBP3-28740](http://www.novusbio.com/reviews/destination/NBP3-28740)



**NBP3-28740-100ug**

ST2/IL-33R Antibody - IgG2SA, (Research Grade melrilimab Biosimilar) - Low Endotoxin, Azide and BSA Free

| Product Information         |                                                                                                                                       |
|-----------------------------|---------------------------------------------------------------------------------------------------------------------------------------|
| Unit Size                   | 100 ug                                                                                                                                |
| Concentration               | LYOPH mg/ml                                                                                                                           |
| Storage                     | Store at 4C.                                                                                                                          |
| Clonality                   | Monoclonal                                                                                                                            |
| Preservative                | No Preservative                                                                                                                       |
| Reconstitution Instructions | Reconstitute with sterile, distilled water to a final concentration of 1 mg/ml. Gently shake to solubilize completely. Do not vortex. |
| Isotype                     | IgG2                                                                                                                                  |
| Purity                      | Protein A purified                                                                                                                    |
| Buffer                      | Lyophilized from 25mM histidine, 8% sucrose, 0.01% Tween80 (pH6.2)                                                                    |

| Product Description |                                                                                                                                                                                                                                                                                                                                                                                                  |
|---------------------|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Description         | <p>Expressed from CHO. The heavy chain type is hulG2, and the light chain type is hukappa. It has a predicted MW of 145.5 kDa.</p> <p>Also known as 'melrilimab'.</p> <p>IgG2SA is human IgG2 with S378A mutation.</p> <p>Upon receipt, store immediately at -20C or lower for 24 months in a lyophilized state. - 80C for 3 months after reconstitution. Avoid repeated freeze-thaw cycles.</p> |
| Host                | Human                                                                                                                                                                                                                                                                                                                                                                                            |
| Gene ID             | 9173                                                                                                                                                                                                                                                                                                                                                                                             |
| Gene Symbol         | IL1RL1                                                                                                                                                                                                                                                                                                                                                                                           |
| Species             | Human                                                                                                                                                                                                                                                                                                                                                                                            |
| Immunogen           | IL-1RL1 / ST2 / IL-33R                                                                                                                                                                                                                                                                                                                                                                           |
| Endotoxin Note      | < 0.001EU/ug,determined by LAL method.                                                                                                                                                                                                                                                                                                                                                           |

| Product Application Details |                                   |
|-----------------------------|-----------------------------------|
| Applications                | ELISA, Flow Cytometry, Functional |
| Recommended Dilutions       | Flow Cytometry, ELISA, Functional |





### **Novus Biologicals USA**

10730 E. Briarwood Avenue  
Centennial, CO 80112  
USA  
Phone: 303.730.1950  
Toll Free: 1.888.506.6887  
Fax: 303.730.1966  
nb-customerservice@bio-techne.com

### **Bio-Techne Canada**

21 Canmotor Ave  
Toronto, ON M8Z 4E6  
Canada  
Phone: 905.827.6400  
Toll Free: 855.668.8722  
Fax: 905.827.6402  
canada.inquires@bio-techne.com

### **Bio-Techne Ltd**

19 Barton Lane  
Abingdon Science Park  
Abingdon, OX14 3NB, United Kingdom  
Phone: (44) (0) 1235 529449  
Free Phone: 0800 37 34 15  
Fax: (44) (0) 1235 533420  
info.EMEA@bio-techne.com

### **General Contact Information**

[www.novusbio.com](http://www.novusbio.com)  
Technical Support: nb-technical@bio-techne.com  
Orders: nb-customerservice@bio-techne.com  
General: novus@novusbio.com

### **Products Related to NBP3-28740-100ug**

|               |                                                       |
|---------------|-------------------------------------------------------|
| G-102-C       | Goat anti-Human IgG Secondary Antibody [Unconjugated] |
| NB7446        | Goat anti-Human IgG Fc Secondary Antibody             |
| DDXCH02P-100  | Human IgG2 Isotype Control                            |
| NBP2-24706PEP | ST2/IL-33R Antibody Blocking Peptide                  |

### **Limitations**

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit [www.novusbio.com/guarantee](http://www.novusbio.com/guarantee)

Earn gift cards/discounts by submitting a review: [www.novusbio.com/reviews/submit/NBP3-28740](http://www.novusbio.com/reviews/submit/NBP3-28740)

Earn gift cards/discounts by submitting a publication using this product:  
[www.novusbio.com/publications](http://www.novusbio.com/publications)

