

# Product Datasheet

## Dynactin Subunit 2/DCTN2/DCTN-50 Antibody NBP3-29034

Unit Size: 20 ug

Store at 4C. Do not freeze.

[www.novusbio.com](http://www.novusbio.com)



[technical@novusbio.com](mailto:technical@novusbio.com)

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:  
[www.novusbio.com/NBP3-29034](http://www.novusbio.com/NBP3-29034)

Updated 9/9/2025 v.20.1

Earn rewards for product  
reviews and publications.

Submit a publication at [www.novusbio.com/publications](http://www.novusbio.com/publications)

Submit a review at [www.novusbio.com/reviews/destination/NBP3-29034](http://www.novusbio.com/reviews/destination/NBP3-29034)



**NBP3-29034****Dynactin Subunit 2/DCTN2/DCTN-50 Antibody****Product Information**

|                      |                                                                                             |
|----------------------|---------------------------------------------------------------------------------------------|
| <b>Unit Size</b>     | 20 ug                                                                                       |
| <b>Concentration</b> | Please see the vial label for concentration. If unlisted please contact technical services. |
| <b>Storage</b>       | Store at 4C. Do not freeze.                                                                 |
| <b>Clonality</b>     | Polyclonal                                                                                  |
| <b>Preservative</b>  | 0.09% Sodium Azide                                                                          |
| <b>Isotype</b>       | IgG                                                                                         |
| <b>Purity</b>        | Affinity purified                                                                           |
| <b>Buffer</b>        | Tris-buffered Saline 0.1% BSA                                                               |

**Product Description**


|                      |                                                                                                                                                                                                                   |
|----------------------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Description</b>   | Novus Biologicals Rabbit Dynactin Subunit 2/DCTN2/DCTN-50 Antibody (NBP3-29034) is a polyclonal antibody validated for use in IHC, WB and IP. All Novus Biologicals antibodies are covered by our 100% guarantee. |
| <b>Host</b>          | Rabbit                                                                                                                                                                                                            |
| <b>Gene ID</b>       | 10540                                                                                                                                                                                                             |
| <b>Gene Symbol</b>   | DCTN2                                                                                                                                                                                                             |
| <b>Species</b>       | Human                                                                                                                                                                                                             |
| <b>Immunogen</b>     | between 351 and 401                                                                                                                                                                                               |
| <b>Assay Type</b>    | [Assay Type]                                                                                                                                                                                                      |
| <b>Sample Volume</b> | [Sample Volume]                                                                                                                                                                                                   |

**Product Application Details**

|                              |                                                         |
|------------------------------|---------------------------------------------------------|
| <b>Applications</b>          | Western Blot, Immunohistochemistry, Immunoprecipitation |
| <b>Recommended Dilutions</b> | Western Blot, Immunohistochemistry, Immunoprecipitation |

**Images**

Immunohistochemistry:[NBP3-29034]-Detection of human DCTN2 by immunohistochemistry. Sample: FFPE section of human breast carcinoma. Antibody: Affinity purified rabbit anti-DCTN2 used at a dilution of 1: 200. Detection: DAB. Counterstain: Hematoxylin.





### **Novus Biologicals USA**

10730 E. Briarwood Avenue  
Centennial, CO 80112  
USA  
Phone: 303.730.1950  
Toll Free: 1.888.506.6887  
Fax: 303.730.1966  
nb-customerservice@bio-techne.com

### **Bio-Techne Canada**

21 Canmotor Ave  
Toronto, ON M8Z 4E6  
Canada  
Phone: 905.827.6400  
Toll Free: 855.668.8722  
Fax: 905.827.6402  
canada.inquires@bio-techne.com

### **Bio-Techne Ltd**

19 Barton Lane  
Abingdon Science Park  
Abingdon, OX14 3NB, United Kingdom  
Phone: (44) (0) 1235 529449  
Free Phone: 0800 37 34 15  
Fax: (44) (0) 1235 533420  
info.EMEA@bio-techne.com

### **General Contact Information**

[www.novusbio.com](http://www.novusbio.com)  
Technical Support: nb-technical@bio-techne.com  
Orders: nb-customerservice@bio-techne.com  
General: novus@novusbio.com

### **Products Related to NBP3-29034**

---

|            |                                                                |
|------------|----------------------------------------------------------------|
| HAF008     | Goat anti-Rabbit IgG Secondary Antibody [HRP]                  |
| NB7160     | Goat anti-Rabbit IgG (H+L) Secondary Antibody [HRP]            |
| NBP2-24891 | Rabbit IgG Isotype Control                                     |
| NBP2-23012 | Recombinant Human Dynactin Subunit 2/DCTN2/DCTN-50 His Protein |

---

### **Limitations**

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit [www.novusbio.com/guarantee](http://www.novusbio.com/guarantee)

Earn gift cards/discounts by submitting a review: [www.novusbio.com/reviews/submit/NBP3-29034](http://www.novusbio.com/reviews/submit/NBP3-29034)

Earn gift cards/discounts by submitting a publication using this product:  
[www.novusbio.com/publications](http://www.novusbio.com/publications)

