

Product Datasheet

HA-1 Antibody - BSA Free NBP3-29796

Unit Size: 100 ug

Store at 4C. Do not freeze.

www.novusbio.com



technical@novusbio.com

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/NBP3-29796

Updated 3/18/2025 v.20.1

Earn rewards for product
reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/NBP3-29796



NBP3-29796

HA-1 Antibody - BSA Free

Product Information

| | |
|----------------------|---------------------------------------------------------------------------------------------|
| Unit Size | 100 ug |
| Concentration | Please see the vial label for concentration. If unlisted please contact technical services. |
| Storage | Store at 4C. Do not freeze. |
| Clonality | Polyclonal |
| Preservative | 0.09% Sodium Azide |
| Isotype | IgG |
| Purity | Affinity purified |
| Buffer | Tris-citrate/phosphate pH 7 -8. |

Product Description


| | |
|----------------------|---------------------|
| Host | Rabbit |
| Gene ID | 23526 |
| Gene Symbol | ARHGAP45 |
| Species | Human |
| Immunogen | Between 100 and 150 |
| Assay Type | [Assay Type] |
| Sample Volume | [Sample Volume] |

Product Application Details

| | |
|------------------------------|-----------------------------------|
| Applications | Western Blot, Immunoprecipitation |
| Recommended Dilutions | Western Blot, Immunoprecipitation |

Images

Immunoprecipitation:HA-1 AntibodyNBP3-29796] -Detection of human HMHA1/HA-1 by western blot of immunoprecipitates. Samples: Whole cell lysate (0.5 or 1.0 mg per IP reaction; 20% of IP loaded) from Jurkat cells prepared using NETN lysis buffer. Antibodies: Affinity purified rabbit anti-HMHA1/HA-1 antibody NBP3-29796 (lot NBP3-29796-1) used for IP at 6 µg per reaction. HMHA1/HA-1 was also immunoprecipitated by rabbit anti-HMHA1/HA-1.. For blotting immunoprecipitated HMHA1/HA-1, NBP3-29796 was used at 0.4 µg/ml. Detection: Chemiluminescence with an exposure time of 30 seconds.





Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: nb-technical@bio-techne.com
Orders: nb-customerservice@bio-techne.com
General: novus@novusbio.com

Products Related to NBP3-29796

| | |
|---------------|-----------------------------------------------------|
| HAF008 | Goat anti-Rabbit IgG Secondary Antibody [HRP] |
| NB7160 | Goat anti-Rabbit IgG (H+L) Secondary Antibody [HRP] |
| NBP2-24891 | Rabbit IgG Isotype Control |
| NBP3-17053PEP | HA-1 Recombinant Protein Antigen |

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/NBP3-29796

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

