

Product Datasheet

ABCF3 Antibody - BSA Free

NBP3-29881

Unit Size: 100 ug

Store at 4C. Do not freeze.

www.novusbio.com



technical@novusbio.com

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/NBP3-29881

Updated 3/18/2025 v.20.1

Earn rewards for product
reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/NBP3-29881



NBP3-29881

ABCF3 Antibody - BSA Free

Product Information

| | |
|----------------------|---|
| Unit Size | 100 ug |
| Concentration | Please see the vial label for concentration. If unlisted please contact technical services. |
| Storage | Store at 4C. Do not freeze. |
| Clonality | Polyclonal |
| Preservative | 0.09% Sodium Azide |
| Isotype | IgG |
| Purity | Affinity purified |
| Buffer | Tris-citrate/phosphate buffer, pH 7 to 8 |

Product Description


| | |
|----------------------|---------------------|
| Host | Rabbit |
| Gene ID | 55324 |
| Gene Symbol | ABCF3 |
| Species | Human, Mouse |
| Immunogen | between 659 and 709 |
| Assay Type | [Assay Type] |
| Sample Volume | [Sample Volume] |

Product Application Details

| | |
|------------------------------|-----------------------------------|
| Applications | Western Blot, Immunoprecipitation |
| Recommended Dilutions | Western Blot, Immunoprecipitation |

Images

Western Blot:[NBP3-29881]-Detection of human and mouse ABCF3 by western blot and immunoprecipitation. Samples: Whole cell lysate from HeLa , HEK293T , and mouse NIH 3T3 cells. Antibodies: Affinity purified rabbit anti-ABCF3 antibody used for WB at 0. 1 µg/ml and 0. 4 µg/ml and used for IP at 3 µg/mg lysate. ABCF3 was less efficiently immunoprecipitated by rabbit anti-ABCF3 antibody , which recognizes an upstream epitope. Detection: Chemiluminescence with exposure times of 30 seconds and 3 minutes.





Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: nb-technical@bio-techne.com
Orders: nb-customerservice@bio-techne.com
General: novus@novusbio.com

Products Related to NBP3-29881

| | |
|--------------------|---|
| HAF008 | Goat anti-Rabbit IgG Secondary Antibody [HRP] |
| NB7160 | Goat anti-Rabbit IgG (H+L) Secondary Antibody [HRP] |
| NBP2-24891 | Rabbit IgG Isotype Control |
| H00055324-P01-10ug | Recombinant Human ABCF3 GST (N-Term) Protein |

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/NBP3-29881

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

