

# Product Datasheet

## CNOT1 Antibody - BSA Free

### NBP3-29963

Unit Size: 100 ug

Store at 4C. Do not freeze.

[www.novusbio.com](http://www.novusbio.com)



[technical@novusbio.com](mailto:technical@novusbio.com)

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:  
[www.novusbio.com/NBP3-29963](http://www.novusbio.com/NBP3-29963)

Updated 9/9/2025 v.20.1

Earn rewards for product  
reviews and publications.

Submit a publication at [www.novusbio.com/publications](http://www.novusbio.com/publications)

Submit a review at [www.novusbio.com/reviews/destination/NBP3-29963](http://www.novusbio.com/reviews/destination/NBP3-29963)



**NBP3-29963**

CNOT1 Antibody - BSA Free

**Product Information**

|                      |                                                                                             |
|----------------------|---------------------------------------------------------------------------------------------|
| <b>Unit Size</b>     | 100 ug                                                                                      |
| <b>Concentration</b> | Please see the vial label for concentration. If unlisted please contact technical services. |
| <b>Storage</b>       | Store at 4C. Do not freeze.                                                                 |
| <b>Clonality</b>     | Polyclonal                                                                                  |
| <b>Preservative</b>  | 0.09% Sodium Azide                                                                          |
| <b>Isotype</b>       | IgG                                                                                         |
| <b>Purity</b>        | Affinity purified                                                                           |
| <b>Buffer</b>        | Tris-citrate/phosphate pH 7 -8.                                                             |

**Product Description**


|                      |                                                                                                                                                                                              |
|----------------------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Description</b>   | Novus Biologicals Rabbit CNOT1 Antibody - BSA Free (NBP3-29963) is a polyclonal antibody validated for use in WB and IP. All Novus Biologicals antibodies are covered by our 100% guarantee. |
| <b>Host</b>          | Rabbit                                                                                                                                                                                       |
| <b>Gene ID</b>       | 23019                                                                                                                                                                                        |
| <b>Gene Symbol</b>   | CNOT1                                                                                                                                                                                        |
| <b>Species</b>       | Human, Mouse                                                                                                                                                                                 |
| <b>Immunogen</b>     | between 1400 and 1450                                                                                                                                                                        |
| <b>Assay Type</b>    | [Assay Type]                                                                                                                                                                                 |
| <b>Sample Volume</b> | [Sample Volume]                                                                                                                                                                              |

**Product Application Details**

|                              |                                   |
|------------------------------|-----------------------------------|
| <b>Applications</b>          | Western Blot, Immunoprecipitation |
| <b>Recommended Dilutions</b> | Western Blot, Immunoprecipitation |

**Images**

Immunoprecipitation:CNOT1 AntibodyNBP3-29963] -Detection of human CNOT1 by western blot of immunoprecipitates. Samples: Whole cell lysate (1.0 mg per IP reaction; 20% of IP loaded) from HeLa cells prepared using NETN lysis buffer. Antibodies: Affinity purified Rabbit anti-CNOT1 antibody NBP3-29963 (lot NBP3-29963-2) used for IP at 100 µl per reaction. CNOT1 was also immunoprecipitated by a previous lot of this antibody (lot NBP3-29963). For blotting immunoprecipitated CNOT1, NBP3-29963 was used at 1:1000. Detection: Chemiluminescence with an exposure time of 3 minutes.





### **Novus Biologicals USA**

10730 E. Briarwood Avenue  
Centennial, CO 80112  
USA  
Phone: 303.730.1950  
Toll Free: 1.888.506.6887  
Fax: 303.730.1966  
nb-customerservice@bio-techne.com

### **Bio-Techne Canada**

21 Canmotor Ave  
Toronto, ON M8Z 4E6  
Canada  
Phone: 905.827.6400  
Toll Free: 855.668.8722  
Fax: 905.827.6402  
canada.inquires@bio-techne.com

### **Bio-Techne Ltd**

19 Barton Lane  
Abingdon Science Park  
Abingdon, OX14 3NB, United Kingdom  
Phone: (44) (0) 1235 529449  
Free Phone: 0800 37 34 15  
Fax: (44) (0) 1235 533420  
info.EMEA@bio-techne.com

### **General Contact Information**

www.novusbio.com  
Technical Support: nb-technical@bio-techne.com  
Orders: nb-customerservice@bio-techne.com  
General: novus@novusbio.com

### **Products Related to NBP3-29963**

---

|                    |                                                     |
|--------------------|-----------------------------------------------------|
| HAF008             | Goat anti-Rabbit IgG Secondary Antibody [HRP]       |
| NB7160             | Goat anti-Rabbit IgG (H+L) Secondary Antibody [HRP] |
| NBP2-24891         | Rabbit IgG Isotype Control                          |
| H00023019-Q01-10ug | Recombinant Human CNOT1 GST (N-Term) Protein        |

---

### **Limitations**

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit [www.novusbio.com/guarantee](http://www.novusbio.com/guarantee)

Earn gift cards/discounts by submitting a review: [www.novusbio.com/reviews/submit/NBP3-29963](http://www.novusbio.com/reviews/submit/NBP3-29963)

Earn gift cards/discounts by submitting a publication using this product:  
[www.novusbio.com/publications](http://www.novusbio.com/publications)

