

Product Datasheet

BLOS3 Antibody - BSA Free NBP3-30074

Unit Size: 100 ug

Store at 4C. Do not freeze.

www.novusbio.com



technical@novusbio.com

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/NBP3-30074

Updated 9/9/2025 v.20.1

Earn rewards for product
reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/NBP3-30074



NBP3-30074

BLOS3 Antibody - BSA Free


| Product Information | |
|---------------------|---|
| Unit Size | 100 ug |
| Concentration | Please see the vial label for concentration. If unlisted please contact technical services. |
| Storage | Store at 4C. Do not freeze. |
| Clonality | Polyclonal |
| Preservative | 0.09% Sodium Azide |
| Isotype | IgG |
| Purity | Affinity purified |
| Buffer | Tris-citrate/phosphate buffer, pH 7 to 8 |

| Product Description | |
|---------------------|--|
| Description | Novus Biologicals Rabbit BLOS3 Antibody - BSA Free (NBP3-30074) is a polyclonal antibody validated for use in WB and IP. All Novus Biologicals antibodies are covered by our 100% guarantee. |
| Host | Rabbit |
| Gene ID | 388552 |
| Gene Symbol | BLOC1S3 |
| Species | Human |
| Immunogen | Between 1 and 50 |
| Assay Type | [Assay Type] |
| Sample Volume | [Sample Volume] |

| Product Application Details | |
|-----------------------------|-----------------------------------|
| Applications | Western Blot, Immunoprecipitation |
| Recommended Dilutions | Western Blot, Immunoprecipitation |

Images

Immunoprecipitation: BLOS3 Antibody NBP3-30074] -Detection of human BLOC1S3 by western blot of immunoprecipitates. Samples: Whole cell lysate (0.5 or 1.0 mg per IP reaction; 20% of IP loaded) from 293T cells prepared using NETN lysis buffer. Antibodies: Affinity purified rabbit anti-BLOC1S3 antibody NBP3-30074 (lot NBP3-30074-1) used for IP at 6 µg per reaction. BLOC1S3 was also immunoprecipitated by rabbit anti-BLOC1S3 antibody BL17617. For blotting immunoprecipitated BLOC1S3, NBP3-30074 was used at 1 µg/ml. Detection: Chemiluminescence with an exposure time of 30 seconds.





Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: nb-technical@bio-techne.com
Orders: nb-customerservice@bio-techne.com
General: novus@novusbio.com

Products Related to NBP3-30074

| | |
|-------------|---|
| NBP2-33376H | Blue Marker Antibody (6F4-F6) [HRP] |
| HAF008 | Goat anti-Rabbit IgG Secondary Antibody [HRP] |
| NB7160 | Goat anti-Rabbit IgG (H+L) Secondary Antibody [HRP] |
| NBP2-24891 | Rabbit IgG Isotype Control |

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/NBP3-30074

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

