

Product Datasheet

C2orf49 Antibody - BSA Free

NBP3-30211

Unit Size: 100 ug

Store at 4C. Do not freeze.

www.novusbio.com



technical@novusbio.com

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/NBP3-30211

Updated 9/9/2025 v.20.1

Earn rewards for product
reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/NBP3-30211



NBP3-30211

C2orf49 Antibody - BSA Free


| Product Information | |
|---------------------|---------------------------------------------------------------------------------------------|
| Unit Size | 100 ug |
| Concentration | Please see the vial label for concentration. If unlisted please contact technical services. |
| Storage | Store at 4C. Do not freeze. |
| Clonality | Polyclonal |
| Preservative | 0.09% Sodium Azide |
| Isotype | IgG |
| Purity | Affinity purified |
| Buffer | Tris-citrate/phosphate pH 7 -8. |

| Product Description | |
|---------------------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Description | Novus Biologicals Rabbit C2orf49 Antibody - BSA Free (NBP3-30211) is a polyclonal antibody validated for use in IP. All Novus Biologicals antibodies are covered by our 100% guarantee. |
| Host | Rabbit |
| Gene ID | 79074 |
| Gene Symbol | C2orf49 |
| Species | Human |
| Immunogen | Between 182 and 232 |
| Assay Type | [Assay Type] |
| Sample Volume | [Sample Volume] |

| Product Application Details | |
|-----------------------------|---------------------|
| Applications | Immunoprecipitation |
| Recommended Dilutions | Immunoprecipitation |

Images

Immunoprecipitation: C2orf49 Antibody [NBP3-30211] - Detection of human C2orf49 by western blot of immunoprecipitates. Samples: Whole cell lysate (0.5 or 1.0 mg per IP reaction; 20% of IP loaded) from 293T cells prepared using NETN lysis buffer. Antibodies: Affinity purified rabbit anti-C2orf49 antibody NBP3-30211 (lot NBP3-30211-1) used for IP at 6 µg per reaction. C2orf49 was also immunoprecipitated by rabbit anti-C2orf49 antibody. For blotting immunoprecipitated C2orf49, NBP3-30211 was used at 0.4 µg/ml. Detection: Chemiluminescence with an exposure time of 3 minutes.





Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: nb-technical@bio-techne.com
Orders: nb-customerservice@bio-techne.com
General: novus@novusbio.com

Products Related to NBP3-30211

| | |
|---------------|-----------------------------------------------------|
| HAF008 | Goat anti-Rabbit IgG Secondary Antibody [HRP] |
| NB7160 | Goat anti-Rabbit IgG (H+L) Secondary Antibody [HRP] |
| NBP2-24891 | Rabbit IgG Isotype Control |
| NBP2-14404PEP | C2orf49 Recombinant Protein Antigen |

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/NBP3-30211

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

