

Product Datasheet

IMP2/IGF2BP2 Antibody (1E3.01E5) NBP3-30645

Unit Size: 100 ul

Store at 4C. Do not freeze.

www.novusbio.com



technical@novusbio.com

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/NBP3-30645

Updated 9/9/2025 v.20.1

Earn rewards for product
reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/NBP3-30645



NBP3-30645

IMP2/IGF2BP2 Antibody (1E3.01E5)

Product Information

Unit Size	100 ul
Concentration	Please see the vial label for concentration. If unlisted please contact technical services.
Storage	Store at 4C. Do not freeze.
Clonality	Monoclonal
Clone	1E3.01E5
Preservative	0.09% Sodium Azide
Isotype	IgG
Purity	Affinity purified
Buffer	Tris-buffered Saline 0.1% BSA

Product Description


Description	Novus Biologicals Mouse IMP2/IGF2BP2 Antibody (1E3.01E5) (NBP3-30645) is a monoclonal antibody validated for use in WB and IP. All Novus Biologicals antibodies are covered by our 100% guarantee.
Host	Mouse
Gene ID	10644
Gene Symbol	IGF2BP2
Species	Human, Mouse
Immunogen	Full length recombinant human protein IGF mRNA-binding protein 2
Assay Type	[Assay Type]
Sample Volume	[Sample Volume]

Product Application Details

Applications	Western Blot, Immunoprecipitation
Recommended Dilutions	Western Blot, Immunoprecipitation

Images

Immunoprecipitation:[NBP3-30645]-Detection of human IGF2BP2 by western blot of immunoprecipitates. Samples: Whole cell lysate from Hep-G2 cells prepared using NETN lysis buffer. Antibodies: Mouse anti-IGF2BP2 monoclonal antibody [1E3. 01E5] used for IP at 8 µl per reaction. IGF2BP2 was also immunoprecipitated by a previous lot of this antibody and mouse anti-IGF2BP2 antibody. For blotting immunoprecipitated IGF2BP2, was used at 1: 1000. Detection: Chemiluminescence with an exposure time of 10 seconds.



**Novus Biologicals USA**

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: nb-technical@bio-techne.com
Orders: nb-customerservice@bio-techne.com
General: novus@novusbio.com

Products Related to NBP3-30645

HAF007	Goat anti-Mouse IgG Secondary Antibody [HRP]
NB720-B	Rabbit anti-Mouse IgG (H+L) Secondary Antibody [Biotin]
NBP1-97019-5mg	Mouse IgG Isotype Control
H00010644-Q01-10ug	Recombinant Human IMP2/IGF2BP2 GST (N-Term) Protein

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/NBP3-30645

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

