

# Product Datasheet

## Human A33/GPA33 - Ready-To-Use ELISA Kit (Colorimetric) NBP3-30759

Unit Size: 1 Kit

Storage of components varies. See protocol for specific instructions.

[www.novusbio.com](http://www.novusbio.com)



[technical@novusbio.com](mailto:technical@novusbio.com)

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:  
[www.novusbio.com/NBP3-30759](http://www.novusbio.com/NBP3-30759)

Updated 10/7/2024 v.20.1

Earn rewards for product  
reviews and publications.

Submit a publication at [www.novusbio.com/publications](http://www.novusbio.com/publications)

Submit a review at [www.novusbio.com/reviews/destination/NBP3-30759](http://www.novusbio.com/reviews/destination/NBP3-30759)




**NBP3-30759****Human A33/GPA33 - Ready-To-Use ELISA Kit (Colorimetric)**

| Product Information         |                                                                                                                                                                                                                 |
|-----------------------------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Unit Size                   | 1 Kit                                                                                                                                                                                                           |
| Concentration               | Concentration is not relevant for this product. Please see the protocols for proper use of this product.                                                                                                        |
| Storage                     | Storage of components varies. See protocol for specific instructions.                                                                                                                                           |
| Product Description         |                                                                                                                                                                                                                 |
| Description                 | The Ready-To-Use ELISA kit offers pre-diluted detection reagents and a shorter experimental time.<br>Assay Length: 3 hours                                                                                      |
| Gene ID                     | 10223                                                                                                                                                                                                           |
| Gene Symbol                 | GPA33                                                                                                                                                                                                           |
| Species                     | Human                                                                                                                                                                                                           |
| Kit Components              | Pre-coated 96T strip plate, Plate sealer for 96 wells, Standard, Standard Diluent, Detection Solution A, Detection Solution B, TMB Substrate, Stop Solution, Wash Buffer (30 x concentrate), Instruction manual |
| Standard Curve Range        | 0.625 - 40 ng/mL (example only; lot dependent)                                                                                                                                                                  |
| Sensitivity                 | 0.238 ng/mL (example only; lot dependent)                                                                                                                                                                       |
| Inter Assay                 | %CV < 12 (example only; lot dependent)                                                                                                                                                                          |
| Intra Assay                 | %CV < 10 (example only; lot dependent)                                                                                                                                                                          |
| Assay Type                  | Sandwich ELISA                                                                                                                                                                                                  |
| Suitable Sample Type        | Tissue Homogenates, Cell Lysates and other Biological Fluids                                                                                                                                                    |
| Sample Volume               | 100 $\mu$ L                                                                                                                                                                                                     |
| Product Application Details |                                                                                                                                                                                                                 |
| Applications                | ELISA                                                                                                                                                                                                           |
| Recommended Dilutions       | ELISA                                                                                                                                                                                                           |

**Images**

ELISA: Human A33/GPA33 - Ready-To-Use ELISA Kit (Colorimetric)  
[NBP3-30759] - Standard Curve Reference





### **Novus Biologicals USA**

10730 E. Briarwood Avenue  
Centennial, CO 80112  
USA  
Phone: 303.730.1950  
Toll Free: 1.888.506.6887  
Fax: 303.730.1966  
nb-customerservice@bio-techne.com

### **Bio-Techne Canada**

21 Canmotor Ave  
Toronto, ON M8Z 4E6  
Canada  
Phone: 905.827.6400  
Toll Free: 855.668.8722  
Fax: 905.827.6402  
canada.inquires@bio-techne.com

### **Bio-Techne Ltd**

19 Barton Lane  
Abingdon Science Park  
Abingdon, OX14 3NB, United Kingdom  
Phone: (44) (0) 1235 529449  
Free Phone: 0800 37 34 15  
Fax: (44) (0) 1235 533420  
info.EMEA@bio-techne.com

### **General Contact Information**

[www.novusbio.com](http://www.novusbio.com)  
Technical Support: nb-technical@bio-techne.com  
Orders: nb-customerservice@bio-techne.com  
General: novus@novusbio.com

### **Products Related to NBP3-30759**

---

|                    |                                                  |
|--------------------|--------------------------------------------------|
| H00010223-P01-10ug | Recombinant Human A33/GPA33 GST (N-Term) Protein |
| 210-TA-005         | TNF-alpha [Unconjugated]                         |
| AF3080             | A33/GPA33 Antibody [Unconjugated]                |
| 314-BP-010         | BMP-4 [Unconjugated]                             |

---

### **Limitations**

This product is for research use only and is not approved for use in humans or in clinical diagnosis. ELISA Kits are guaranteed for 6 months from date of receipt.

For more information on our 100% guarantee, please visit [www.novusbio.com/guarantee](http://www.novusbio.com/guarantee)

Earn gift cards/discounts by submitting a review: [www.novusbio.com/reviews/submit/NBP3-30759](http://www.novusbio.com/reviews/submit/NBP3-30759)

Earn gift cards/discounts by submitting a publication using this product:  
[www.novusbio.com/publications](http://www.novusbio.com/publications)

