

# Product Datasheet

## Human Glucosylceramidase/GBA - Ready-To-Use ELISA Kit (Colorimetric) NBP3-31768

Unit Size: 1 Kit

Storage of components varies. See protocol for specific instructions.

[www.novusbio.com](http://www.novusbio.com)



[technical@novusbio.com](mailto:technical@novusbio.com)

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:  
[www.novusbio.com/NBP3-31768](http://www.novusbio.com/NBP3-31768)

Updated 6/19/2024 v.20.1

Earn rewards for product  
reviews and publications.

Submit a publication at [www.novusbio.com/publications](http://www.novusbio.com/publications)

Submit a review at [www.novusbio.com/reviews/destination/NBP3-31768](http://www.novusbio.com/reviews/destination/NBP3-31768)



**NBP3-31768****Human Glucosylceramidase/GBA - Ready-To-Use ELISA Kit (Colorimetric)**

| <b>Product Information</b>         |                                                                                                                                                                                                                 |
|------------------------------------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Unit Size</b>                   | 1 Kit                                                                                                                                                                                                           |
| <b>Concentration</b>               | Concentration is not relevant for this product. Please see the protocols for proper use of this product.                                                                                                        |
| <b>Storage</b>                     | Storage of components varies. See protocol for specific instructions.                                                                                                                                           |
| <b>Product Description</b>         |                                                                                                                                                                                                                 |
| <b>Description</b>                 | The Ready-To-Use ELISA kit offers pre-diluted detection reagents and a shorter experimental time.<br>Assay Length: 3 hours                                                                                      |
| <b>Gene ID</b>                     | 2629                                                                                                                                                                                                            |
| <b>Gene Symbol</b>                 | GBA1                                                                                                                                                                                                            |
| <b>Species</b>                     | Human                                                                                                                                                                                                           |
| <b>Kit Components</b>              | Pre-coated 96T strip plate, Plate sealer for 96 wells, Standard, Standard Diluent, Detection Solution A, Detection Solution B, TMB Substrate, Stop Solution, Wash Buffer (30 x concentrate), Instruction manual |
| <b>Standard Curve Range</b>        | 1.56 - 100 ng/mL (example only; lot dependent)                                                                                                                                                                  |
| <b>Sensitivity</b>                 | 0.59 ng/mL (example only; lot dependent)                                                                                                                                                                        |
| <b>Inter Assay</b>                 | %CV < 12 (example only; lot dependent)                                                                                                                                                                          |
| <b>Intra Assay</b>                 | %CV < 10 (example only; lot dependent)                                                                                                                                                                          |
| <b>Assay Type</b>                  | Sandwich ELISA                                                                                                                                                                                                  |
| <b>Suitable Sample Type</b>        | Serum, plasma, tissue homogenates and other biological fluids                                                                                                                                                   |
| <b>Sample Volume</b>               | 100 uL                                                                                                                                                                                                          |
| <b>Product Application Details</b> |                                                                                                                                                                                                                 |
| <b>Applications</b>                | ELISA                                                                                                                                                                                                           |
| <b>Recommended Dilutions</b>       | ELISA                                                                                                                                                                                                           |





### **Novus Biologicals USA**

10730 E. Briarwood Avenue  
Centennial, CO 80112  
USA  
Phone: 303.730.1950  
Toll Free: 1.888.506.6887  
Fax: 303.730.1966  
nb-customerservice@bio-techne.com

### **Bio-Techne Canada**

21 Canmotor Ave  
Toronto, ON M8Z 4E6  
Canada  
Phone: 905.827.6400  
Toll Free: 855.668.8722  
Fax: 905.827.6402  
canada.inquires@bio-techne.com

### **Bio-Techne Ltd**

19 Barton Lane  
Abingdon Science Park  
Abingdon, OX14 3NB, United Kingdom  
Phone: (44) (0) 1235 529449  
Free Phone: 0800 37 34 15  
Fax: (44) (0) 1235 533420  
info.EMEA@bio-techne.com

### **General Contact Information**

[www.novusbio.com](http://www.novusbio.com)  
Technical Support: nb-technical@bio-techne.com  
Orders: nb-customerservice@bio-techne.com  
General: novus@novusbio.com

### **Products Related to NBP3-31768**

---

|                   |                                                      |
|-------------------|------------------------------------------------------|
| NBP2-52148-0.05mg | Recombinant Human Glucosylceramidase/GBA His Protein |
| M6000B-1          | IL-6 [HRP]                                           |
| H00002629-M01     | Glucosylceramidase/GBA Antibody (2E2)                |
| NB300-268         | LRRK2 Antibody - BSA Free                            |

---

### **Limitations**

This product is for research use only and is not approved for use in humans or in clinical diagnosis. ELISA Kits are guaranteed for 6 months from date of receipt.

For more information on our 100% guarantee, please visit [www.novusbio.com/guarantee](http://www.novusbio.com/guarantee)

Earn gift cards/discounts by submitting a review: [www.novusbio.com/reviews/submit/NBP3-31768](http://www.novusbio.com/reviews/submit/NBP3-31768)

Earn gift cards/discounts by submitting a publication using this product:  
[www.novusbio.com/publications](http://www.novusbio.com/publications)

