

Product Datasheet

Hexokinase 2 Antibody (5Z8C3)

NBP3-33349-100ul

Unit Size: 100 ul

Store at -20C. Avoid freeze-thaw cycles.

www.novusbio.com



technical@novusbio.com

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/NBP3-33349

Updated 9/9/2025 v.20.1

Earn rewards for product
reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/NBP3-33349



NBP3-33349-100ul

Hexokinase 2 Antibody (5Z8C3)

Product Information

Unit Size	100 ul
Concentration	Please see the vial label for concentration. If unlisted please contact technical services.
Storage	Store at -20C. Avoid freeze-thaw cycles.
Clonality	Monoclonal
Clone	5Z8C3
Preservative	0.05% Proclin 300
Isotype	IgG
Purity	Affinity purified
Buffer	PBS (pH 7.3), 50% glycerol, 0.05% BSA
Target Molecular Weight	102 kDa

Product Description


Description	Novus Biologicals Rabbit Hexokinase 2 Antibody (5Z8C3) (NBP3-33349) is a monoclonal antibody validated for use in WB and ELISA. All Novus Biologicals antibodies are covered by our 100% guarantee.
Host	Rabbit
Gene ID	3099
Gene Symbol	HK2
Species	Human, Mouse, Rat
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-120 of human Hexokinase 2 (NP_000180.2). Sequence: MIASHLLAYFFTELNHDQVQKVDQYLYHMRLSDETLLEISKRFKEMEKGLGAT THPTAAVKMLPTFVRSTPDGTEHGEFLALDLGGTNFRVLWVKVTDNGLQKVE MENQIYAIPEDIMR

Product Application Details


Applications	Western Blot, ELISA
Recommended Dilutions	Western Blot 1:5000 - 1:30000, ELISA Recommended starting concentration is 1 ug/mL

Images


Western Blot: Hexokinase 2 Antibody (5Z8C3) [NBP3-33349] - Western blot analysis of various lysates using [KO Validated] Hexokinase 2 Rabbit mAb at 1:25000 dilution incubated overnight at 4C. Secondary antibody: HRP-conjugated Goat anti-Rabbit IgG (H+L) at 1:10000 dilution.
Lysates/proteins: 25 ug per lane.
Blocking buffer: 3% nonfat dry milk in TBST.
Detection: ECL Basic Kit.
Exposure time: 30s.




Western Blot: Hexokinase 2 Antibody (5Z8C3) [NBP3-33349] - Western blot analysis of lysates from Mouse skeletal muscle using [KO Validated] Hexokinase 2 Rabbit mAb at 1:25000 dilution incubated overnight at 4C. Secondary antibody: HRP-conjugated Goat anti-Rabbit IgG (H+L) at 1:10000 dilution.
Lysates/proteins: 25 ug per lane.
Blocking buffer: 3% nonfat dry milk in TBST.
Detection: ECL Basic Kit.
Exposure time: 10s.




Western Blot: Hexokinase 2 Antibody (5Z8C3) [NBP3-33349] - Western blot analysis of various lysates using [KO Validated] Hexokinase 2 Rabbit mAb at 1:25000 dilution incubated overnight at 4C. Secondary antibody: HRP-conjugated Goat anti-Rabbit IgG (H+L) at 1:10000 dilution.
Lysates/proteins: 25 ug per lane.
Blocking buffer: 3% nonfat dry milk in TBST.
Detection: ECL Basic Kit.
Exposure time: 30s.



Western Blot: Hexokinase 2 Antibody (5Z8C3) [NBP3-33349] - Western blot analysis of lysates from Mouse skeletal muscle using [KO Validated] Hexokinase 2 Rabbit mAb at 1:25000 dilution incubated overnight at 4C. Secondary antibody: HRP-conjugated Goat anti-Rabbit IgG (H+L) at 1:10000 dilution.
Lysates/proteins: 25 ug per lane.
Blocking buffer: 3% nonfat dry milk in TBST.
Detection: ECL Basic Kit.
Exposure time: 10s.



Western Blot: Hexokinase 2 Antibody (5Z8C3) [Hexokinase 2] - Western blot analysis of various lysates using [KO Validated] Hexokinase 2 Rabbit mAb at 1:25000 dilution incubated overnight at 4C. Secondary antibody: HRP-conjugated Goat anti-Rabbit IgG (H+L) at 1:10000 dilution.
Lysates/proteins: 25 ug per lane.
Blocking buffer: 3% nonfat dry milk in TBST.
Detection: ECL Basic Kit .
Exposure time: 30s.




Western Blot: Hexokinase 2 Antibody (5Z8C3) [Hexokinase 2] - Western blot analysis of lysates from Mouse skeletal muscle using [KO Validated] Hexokinase 2 Rabbit mAb at 1:25000 dilution incubated overnight at 4C. Secondary antibody: HRP-conjugated Goat anti-Rabbit IgG (H+L) at 1:10000 dilution.

Lysates/proteins: 25 ug per lane.

Blocking buffer: 3% nonfat dry milk in TBST.

Detection: ECL Basic Kit .

Exposure time: 10s.





Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: nb-technical@bio-techne.com
Orders: nb-customerservice@bio-techne.com
General: novus@novusbio.com

Products Related to NBP3-33349-100ul

HAF008	Goat anti-Rabbit IgG Secondary Antibody [HRP]
NB7160	Goat anti-Rabbit IgG (H+L) Secondary Antibody [HRP]
NBP2-24891	Rabbit IgG Isotype Control
NBP3-21294PEP	Hexokinase 2 Recombinant Protein Antigen

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/NBP3-33349

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

