

# Product Datasheet

## Pea3 Antibody (10I1V4) NBP3-33399-100ul

Unit Size: 100 ul

Store at -20C. Avoid freeze-thaw cycles.

[www.novusbio.com](http://www.novusbio.com)



[technical@novusbio.com](mailto:technical@novusbio.com)

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:  
[www.novusbio.com/NBP3-33399](http://www.novusbio.com/NBP3-33399)

Updated 10/15/2024 v.20.1

Earn rewards for product  
reviews and publications.

Submit a publication at [www.novusbio.com/publications](http://www.novusbio.com/publications)

Submit a review at [www.novusbio.com/reviews/destination/NBP3-33399](http://www.novusbio.com/reviews/destination/NBP3-33399)



**NBP3-33399-100ul**

Pea3 Antibody (10I1V4)

| Product Information         |                                                                                                                                                                                                                                                 |
|-----------------------------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Unit Size                   | 100 ul                                                                                                                                                                                                                                          |
| Concentration               | Please see the vial label for concentration. If unlisted please contact technical services.                                                                                                                                                     |
| Storage                     | Store at -20C. Avoid freeze-thaw cycles.                                                                                                                                                                                                        |
| Clonality                   | Monoclonal                                                                                                                                                                                                                                      |
| Clone                       | 10I1V4                                                                                                                                                                                                                                          |
| Preservative                | 0.05% Proclin 300                                                                                                                                                                                                                               |
| Isotype                     | IgG                                                                                                                                                                                                                                             |
| Purity                      | Affinity purified                                                                                                                                                                                                                               |
| Buffer                      | PBS (pH 7.3), 50% glycerol, 0.05% BSA                                                                                                                                                                                                           |
| Target Molecular Weight     | 54 kDa                                                                                                                                                                                                                                          |
| Product Description         |                                                                                                                                                                                                                                                 |
| Host                        | Rabbit                                                                                                                                                                                                                                          |
| Gene ID                     | 2118                                                                                                                                                                                                                                            |
| Gene Symbol                 | ETV4                                                                                                                                                                                                                                            |
| Species                     | Human                                                                                                                                                                                                                                           |
| Immunogen                   | <p>A synthetic peptide corresponding to a sequence within amino acids 384-484 of human Pea3 (NP_001977.1).</p> <p>Sequence:<br/>QKNRPAMNYDKLSRSLRYYYEKGIMQKVAGERYVYKFVCEPEALFSLAFPDNQ<br/>RPALKAEFDRPVSEEDTVPLSHLDESPAYLPELAGPAQPFPGPKGGYSY</p> |
| Product Application Details |                                                                                                                                                                                                                                                 |
| Applications                | Western Blot, ELISA                                                                                                                                                                                                                             |
| Recommended Dilutions       | Western Blot 1:2000 - 1:9000, ELISA Recommended starting concentration is 1 ug/mL                                                                                                                                                               |





### **Novus Biologicals USA**

10730 E. Briarwood Avenue  
Centennial, CO 80112  
USA  
Phone: 303.730.1950  
Toll Free: 1.888.506.6887  
Fax: 303.730.1966  
nb-customerservice@bio-techne.com

### **Bio-Techne Canada**

21 Canmotor Ave  
Toronto, ON M8Z 4E6  
Canada  
Phone: 905.827.6400  
Toll Free: 855.668.8722  
Fax: 905.827.6402  
canada.inquires@bio-techne.com

### **Bio-Techne Ltd**

19 Barton Lane  
Abingdon Science Park  
Abingdon, OX14 3NB, United Kingdom  
Phone: (44) (0) 1235 529449  
Free Phone: 0800 37 34 15  
Fax: (44) (0) 1235 533420  
info.EMEA@bio-techne.com

### **General Contact Information**

www.novusbio.com  
Technical Support: nb-technical@bio-techne.com  
Orders: nb-customerservice@bio-techne.com  
General: novus@novusbio.com

### **Products Related to NBP3-33399-100ul**

---

|                    |                                                     |
|--------------------|-----------------------------------------------------|
| HAF008             | Goat anti-Rabbit IgG Secondary Antibody [HRP]       |
| NB7160             | Goat anti-Rabbit IgG (H+L) Secondary Antibody [HRP] |
| NBP2-24891         | Rabbit IgG Isotype Control                          |
| H00002118-P01-10ug | Recombinant Human Pea3 GST (N-Term) Protein         |

---

### **Limitations**

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit [www.novusbio.com/guarantee](http://www.novusbio.com/guarantee)

Earn gift cards/discounts by submitting a review: [www.novusbio.com/reviews/submit/NBP3-33399](http://www.novusbio.com/reviews/submit/NBP3-33399)

Earn gift cards/discounts by submitting a publication using this product:  
[www.novusbio.com/publications](http://www.novusbio.com/publications)

