

# Product Datasheet

## RPL32 Antibody NBP3-35299-100ul

Unit Size: 100 ul

Store at -20C. Avoid freeze-thaw cycles.

[www.novusbio.com](http://www.novusbio.com)



[technical@novusbio.com](mailto:technical@novusbio.com)

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:  
[www.novusbio.com/NBP3-35299](http://www.novusbio.com/NBP3-35299)

Updated 9/19/2024 v.20.1

Earn rewards for product  
reviews and publications.

Submit a publication at [www.novusbio.com/publications](http://www.novusbio.com/publications)

Submit a review at [www.novusbio.com/reviews/destination/NBP3-35299](http://www.novusbio.com/reviews/destination/NBP3-35299)



**NBP3-35299-100ul****RPL32 Antibody**

| <b>Product Information</b>         |                                                                                                                                                                                                                                                                                                     |
|------------------------------------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Unit Size</b>                   | 100 ul                                                                                                                                                                                                                                                                                              |
| <b>Concentration</b>               | Please see the vial label for concentration. If unlisted please contact technical services.                                                                                                                                                                                                         |
| <b>Storage</b>                     | Store at -20C. Avoid freeze-thaw cycles.                                                                                                                                                                                                                                                            |
| <b>Clonality</b>                   | Polyclonal                                                                                                                                                                                                                                                                                          |
| <b>Preservative</b>                | 0.01% Thimerosal                                                                                                                                                                                                                                                                                    |
| <b>Isotype</b>                     | IgG                                                                                                                                                                                                                                                                                                 |
| <b>Purity</b>                      | Affinity purified                                                                                                                                                                                                                                                                                   |
| <b>Buffer</b>                      | PBS (pH 7.3), 50% glycerol                                                                                                                                                                                                                                                                          |
| <b>Target Molecular Weight</b>     | 16 kDa                                                                                                                                                                                                                                                                                              |
| <b>Product Description</b>         |                                                                                                                                                                                                                                                                                                     |
| <b>Host</b>                        | Rabbit                                                                                                                                                                                                                                                                                              |
| <b>Gene ID</b>                     | 6161                                                                                                                                                                                                                                                                                                |
| <b>Gene Symbol</b>                 | RPL32                                                                                                                                                                                                                                                                                               |
| <b>Species</b>                     | Human, Mouse, Rat                                                                                                                                                                                                                                                                                   |
| <b>Immunogen</b>                   | <p>Recombinant fusion protein containing a sequence corresponding to amino acids 1-135 of human RPL32 (NP_000985.1).</p> <p>Sequence:<br/> MAALRPLVKPKIVKKRTKKFIRHQSDRYVKIKRNWRKPRGIDNRVRRRFKQGQIL<br/> MPNIGYGSNKKTGHMLPSGFRKFLVHNVKELEVLLMCNKSYPACIAHNVSSKN<br/> RKAIVERAAQLAIRVTNPNARLRSEENE</p> |
| <b>Product Application Details</b> |                                                                                                                                                                                                                                                                                                     |
| <b>Applications</b>                | Western Blot, ELISA, Immunocytochemistry/ Immunofluorescence, Immunohistochemistry, Immunohistochemistry-Paraffin                                                                                                                                                                                   |
| <b>Recommended Dilutions</b>       | Western Blot 1:500 - 1:2000, ELISA, Immunohistochemistry, Immunocytochemistry/ Immunofluorescence 1:50 - 1:100, Immunohistochemistry-Paraffin 1:50 - 1:100                                                                                                                                          |





### **Novus Biologicals USA**

10730 E. Briarwood Avenue  
Centennial, CO 80112  
USA  
Phone: 303.730.1950  
Toll Free: 1.888.506.6887  
Fax: 303.730.1966  
nb-customerservice@bio-techne.com

### **Bio-Techne Canada**

21 Canmotor Ave  
Toronto, ON M8Z 4E6  
Canada  
Phone: 905.827.6400  
Toll Free: 855.668.8722  
Fax: 905.827.6402  
canada.inquires@bio-techne.com

### **Bio-Techne Ltd**

19 Barton Lane  
Abingdon Science Park  
Abingdon, OX14 3NB, United Kingdom  
Phone: (44) (0) 1235 529449  
Free Phone: 0800 37 34 15  
Fax: (44) (0) 1235 533420  
info.EMEA@bio-techne.com

### **General Contact Information**

[www.novusbio.com](http://www.novusbio.com)  
Technical Support: [nb-technical@bio-techne.com](mailto:nb-technical@bio-techne.com)  
Orders: [nb-customerservice@bio-techne.com](mailto:nb-customerservice@bio-techne.com)  
General: [novus@novusbio.com](mailto:novus@novusbio.com)

### **Products Related to NBP3-35299-100ul**

---

|               |                                                     |
|---------------|-----------------------------------------------------|
| HAF008        | Goat anti-Rabbit IgG Secondary Antibody [HRP]       |
| NB7160        | Goat anti-Rabbit IgG (H+L) Secondary Antibody [HRP] |
| NBP2-24891    | Rabbit IgG Isotype Control                          |
| NBP2-30682PEP | RPL32 Recombinant Protein Antigen                   |

---

### **Limitations**

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit [www.novusbio.com/guarantee](http://www.novusbio.com/guarantee)

Earn gift cards/discounts by submitting a review: [www.novusbio.com/reviews/submit/NBP3-35299](http://www.novusbio.com/reviews/submit/NBP3-35299)

Earn gift cards/discounts by submitting a publication using this product:  
[www.novusbio.com/publications](http://www.novusbio.com/publications)

