

Product Datasheet

TRIM27 Antibody - BSA Free NBP3-35634-100ul

Unit Size: 100 ul

Store at -20C. Avoid freeze-thaw cycles.

www.novusbio.com



technical@novusbio.com

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/NBP3-35634

Updated 9/9/2025 v.20.1

Earn rewards for product
reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/NBP3-35634



NBP3-35634-100ul

TRIM27 Antibody - BSA Free

Product Information	
Unit Size	100 ul
Concentration	Please see the vial label for concentration. If unlisted please contact technical services.
Storage	Store at -20C. Avoid freeze-thaw cycles.
Clonality	Polyclonal
Preservative	0.01% Thimerosal
Isotype	IgG
Purity	Affinity purified
Buffer	PBS (pH 7.3), 50% glycerol
Target Molecular Weight	58 kDa


Product Description	
Description	Novus Biologicals Rabbit TRIM27 Antibody - BSA Free (NBP3-35634) is a polyclonal antibody validated for use in IHC and ELISA. All Novus Biologicals antibodies are covered by our 100% guarantee.
Host	Rabbit
Gene ID	5987
Gene Symbol	TRIM27
Species	Human, Mouse, Rat
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-200 of human TRIM27 (NP_006501.1). Sequence: MASGSVAECLQQETTCPVCLQYFAEPMMLDCGHNICCACLARCWGTAETNVS CPQCRETFPQRHMRPNRHLANVTQLVKQLRTERPSGPGGEMGVCEKHREPL KLYCEEDQMPICVVCDRSREHRGHSVLPLEEAVEGFKEQIQNQLDHLKRVKDL KKRRRAQGEQARAELLSLTQMEREKIVWEFEQLYHSLKEHEYRL

Product Application Details	
Applications	Immunohistochemistry-Paraffin, ELISA, Immunohistochemistry
Recommended Dilutions	ELISA, Immunohistochemistry, Immunohistochemistry-Paraffin 1:50 - 1:200




Images


Immunohistochemistry: TRIM27 Antibody [NBP3-35634] -
Immunohistochemistry analysis of paraffin-embedded Human liver using
TRIM27 Rabbit pAb at dilution of 1:100 (40x lens). Microwave antigen
retrieval performed with 0.01M PBS Buffer (pH 7.2) prior to IHC staining.



Immunohistochemistry: TRIM27 Antibody [NBP3-35634] -
Immunohistochemistry analysis of paraffin-embedded Mouse kidney
using TRIM27 Rabbit pAb at dilution of 1:100 (40x lens). Microwave
antigen retrieval performed with 0.01M PBS Buffer (pH 7.2) prior to IHC
staining.



Immunohistochemistry: TRIM27 Antibody [NBP3-35634] -
Immunohistochemistry analysis of paraffin-embedded Mouse kidney
using TRIM27 Rabbit pAb at dilution of 1:100 (40x lens). Perform
microwave antigen retrieval with 10 mM PBS buffer pH 7.2 before
commencing with IHC staining protocol.





Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: nb-technical@bio-techne.com
Orders: nb-customerservice@bio-techne.com
General: novus@novusbio.com

Products Related to NBP3-35634-100ul

HAF008	Goat anti-Rabbit IgG Secondary Antibody [HRP]
NB7160	Goat anti-Rabbit IgG (H+L) Secondary Antibody [HRP]
NBP2-24891	Rabbit IgG Isotype Control
NBP2-54993PEP	TRIM27 Recombinant Protein Antigen

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/NBP3-35634

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

