

Product Datasheet

WDR77 Antibody - BSA Free

NBP3-37891-100ul

Unit Size: 100 ul

Store at -20C. Avoid freeze-thaw cycles.

www.novusbio.com



technical@novusbio.com

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/NBP3-37891

Updated 9/9/2025 v.20.1

Earn rewards for product
reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/NBP3-37891



NBP3-37891-100ul

WDR77 Antibody - BSA Free

Product Information

Unit Size	100 ul
Concentration	Please see the vial label for concentration. If unlisted please contact technical services.
Storage	Store at -20C. Avoid freeze-thaw cycles.
Clonality	Polyclonal
Preservative	0.02% Sodium Azide
Isotype	IgG
Purity	Affinity purified
Buffer	PBS (pH 7.3), 50% glycerol
Target Molecular Weight	43 kDa

Product Description


Description	Novus Biologicals Rabbit WDR77 Antibody - BSA Free (NBP3-37891) is a polyclonal antibody validated for use in IHC, WB, ELISA, ICC/IF and IP. All Novus Biologicals antibodies are covered by our 100% guarantee.
Host	Rabbit
Gene ID	79084
Species	Human
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 290-390 of mouse WDR77 (NP_035770.2). Sequence: PWPLSSFVPSQKTYQGNYGFHLGLQSGTAKSVMCTYSPPLNKLFCQLAKTC PVQLWVSATPPAGSRVRAMAIYKKSQHMTEVVRRCPPHHERCSDGDGLAPPQ HLIRVEGNLYPEYLEDQRQTFRHSVVVPYEPPEAGSEYTTIHYKYMCSNCSMG MNRRPILTIITLEDSSGNLLGRDSFEVRVCACPGRDRTTEENFRKKEVL

Product Application Details


Applications	Western Blot, Immunohistochemistry-Paraffin, ELISA, Immunocytochemistry/Immunofluorescence, Immunohistochemistry, Immunoprecipitation
Recommended Dilutions	Western Blot 1:500 - 1:1000, ELISA, Immunohistochemistry, Immunocytochemistry/ Immunofluorescence 1:50 - 1:200, Immunoprecipitation 0.5ug - 4ug antibody for 200ug - 400ug extracts of whole cells, Immunohistochemistry-Paraffin 1:50 - 1:200

Images


Immunocytochemistry/ Immunofluorescence: WDR77 Antibody [NBP3-37891] - Immunofluorescence analysis of HeLa cells using WDR77 Rabbit pAb. Secondary antibody: Cy3-conjugated Goat anti-Rabbit IgG (H+L) at 1:500 dilution.




Immunoprecipitation: WDR77 Antibody [NBP3-37891] -
 Immunoprecipitation analysis of 200 ug extracts of 293T cells using 3 ug WDR77 antibody. Western blot was performed from the immunoprecipitate using WDR77 antibody at a dilution of 1:1000.




Immunohistochemistry: WDR77 Antibody [NBP3-37891] -
 Immunohistochemistry analysis of paraffin-embedded Human esophageal using WDR77 Rabbit pAb at dilution of 1:100 (40x lens). Microwave antigen retrieval performed with 0.01M PBS Buffer (pH 7.2) prior to IHC staining.




Immunohistochemistry: WDR77 Antibody [NBP3-37891] -
 Immunohistochemistry analysis of paraffin-embedded Human colon carcinoma using WDR77 Rabbit pAb at dilution of 1:100 (40x lens). Microwave antigen retrieval performed with 0.01M PBS Buffer (pH 7.2) prior to IHC staining.




Immunohistochemistry: WDR77 Antibody [NBP3-37891] -
 Immunohistochemistry analysis of paraffin-embedded Human colon using WDR77 Rabbit pAb at dilution of 1:100 (40x lens). Microwave antigen retrieval performed with 0.01M PBS Buffer (pH 7.2) prior to IHC staining.



Immunohistochemistry: WDR77 Antibody [NBP3-37891] - Immunohistochemistry analysis of paraffin-embedded Human esophageal cancer using WDR77 Rabbit pAb at dilution of 1:100 (40x lens). Microwave antigen retrieval performed with 0.01M PBS Buffer (pH 7.2) prior to IHC staining.



Western Blot: WDR77 Antibody [NBP3-37891] - Western blot analysis of lysates from 293T cells, using WDR77 Rabbit pAb at 1:1000 dilution. Secondary antibody: HRP-conjugated Goat anti-Rabbit IgG (H+L) at 1:10000 dilution.
 Lysates/proteins: 25ug per lane.
 Blocking buffer: 3% nonfat dry milk in TBST.
 Detection: ECL Basic Kit.
 Exposure time: 5s.





Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: nb-technical@bio-techne.com
Orders: nb-customerservice@bio-techne.com
General: novus@novusbio.com

Products Related to NBP3-37891-100ul

HAF008	Goat anti-Rabbit IgG Secondary Antibody [HRP]
NB7160	Goat anti-Rabbit IgG (H+L) Secondary Antibody [HRP]
NBP2-24891	Rabbit IgG Isotype Control
H00079084-P01-10ug	Recombinant Human WDR77 GST (N-Term) Protein

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/NBP3-37891

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

