

# Product Datasheet

## Staufen Antibody [DyLight 650] NBP3-38076C

Unit Size: 0.1 ml

Store at 4C in the dark.

[www.novusbio.com](http://www.novusbio.com)



[technical@novusbio.com](mailto:technical@novusbio.com)

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:  
[www.novusbio.com/NBP3-38076C](http://www.novusbio.com/NBP3-38076C)

Updated 3/24/2025 v.20.1

Earn rewards for product  
reviews and publications.

Submit a publication at [www.novusbio.com/publications](http://www.novusbio.com/publications)

Submit a review at [www.novusbio.com/reviews/destination/NBP3-38076C](http://www.novusbio.com/reviews/destination/NBP3-38076C)



**NBP3-38076C**

Staufen Antibody [DyLight 650]

| Product Information          |                                                                                                                                                                                                                                                                                                            |
|------------------------------|------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Unit Size</b>             | 0.1 ml                                                                                                                                                                                                                                                                                                     |
| <b>Concentration</b>         | Please see the vial label for concentration. If unlisted please contact technical services.                                                                                                                                                                                                                |
| <b>Storage</b>               | Store at 4C in the dark.                                                                                                                                                                                                                                                                                   |
| <b>Clonality</b>             | Polyclonal                                                                                                                                                                                                                                                                                                 |
| <b>Preservative</b>          | 0.05% Sodium Azide                                                                                                                                                                                                                                                                                         |
| <b>Isotype</b>               | IgG                                                                                                                                                                                                                                                                                                        |
| <b>Conjugate</b>             | DyLight 650                                                                                                                                                                                                                                                                                                |
| <b>Purity</b>                | Affinity purified                                                                                                                                                                                                                                                                                          |
| <b>Buffer</b>                | 50mM Sodium Borate                                                                                                                                                                                                                                                                                         |
| Product Description          |                                                                                                                                                                                                                                                                                                            |
| <b>Host</b>                  | Rabbit                                                                                                                                                                                                                                                                                                     |
| <b>Gene ID</b>               | 6780                                                                                                                                                                                                                                                                                                       |
| <b>Gene Symbol</b>           | STAU1                                                                                                                                                                                                                                                                                                      |
| <b>Species</b>               | Human, Mouse, Rat                                                                                                                                                                                                                                                                                          |
| <b>Immunogen</b>             | Recombinant fusion protein containing a sequence corresponding to amino acids 426-577 of human Staufen (NP_059347.2).<br><br>Sequence:<br>QAVGVSQGHHTKDFTRAAPNPAKATVTAMIARELLYGGTSPTAETILKNNISSG<br>HVPHGPLTRPSEQLDYLSRVQGFQVEYKDFPKNNKNEFVSLINCSSQPPLISH<br>GIGKDVESCHDMAALNILKLLSELDQQSTEMPRTGNGPMSVCGRC |
| <b>Notes</b>                 | DyLight (R) is a trademark of Thermo Fisher Scientific Inc. and its subsidiaries.                                                                                                                                                                                                                          |
| Product Application Details  |                                                                                                                                                                                                                                                                                                            |
| <b>Applications</b>          | Western Blot, ELISA, Immunohistochemistry, Immunohistochemistry-Paraffin                                                                                                                                                                                                                                   |
| <b>Recommended Dilutions</b> | Western Blot, ELISA, Immunohistochemistry, Immunohistochemistry-Paraffin                                                                                                                                                                                                                                   |
| <b>Application Notes</b>     | Optimal dilution of this antibody should be experimentally determined.                                                                                                                                                                                                                                     |





### **Novus Biologicals USA**

10730 E. Briarwood Avenue  
Centennial, CO 80112  
USA  
Phone: 303.730.1950  
Toll Free: 1.888.506.6887  
Fax: 303.730.1966  
nb-customerservice@bio-techne.com

### **Bio-Techne Canada**

21 Canmotor Ave  
Toronto, ON M8Z 4E6  
Canada  
Phone: 905.827.6400  
Toll Free: 855.668.8722  
Fax: 905.827.6402  
canada.inquires@bio-techne.com

### **Bio-Techne Ltd**

19 Barton Lane  
Abingdon Science Park  
Abingdon, OX14 3NB, United Kingdom  
Phone: (44) (0) 1235 529449  
Free Phone: 0800 37 34 15  
Fax: (44) (0) 1235 533420  
info.EMEA@bio-techne.com

### **General Contact Information**

www.novusbio.com  
Technical Support: nb-technical@bio-techne.com  
Orders: nb-customerservice@bio-techne.com  
General: novus@novusbio.com

### **Products Related to NBP3-38076C**

---

|               |                                          |
|---------------|------------------------------------------|
| NBP2-24891C   | Rabbit IgG Isotype Control [DyLight 650] |
| NBP3-17204PEP | Staufen Recombinant Protein Antigen      |
| NBL1-16531    | Staufen Overexpression Lysate            |
| NB300-213     | MAP2 Antibody - BSA Free                 |

---

### **Limitations**

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit [www.novusbio.com/guarantee](http://www.novusbio.com/guarantee)

Earn gift cards/discounts by submitting a review: [www.novusbio.com/reviews/submit/NBP3-38076C](http://www.novusbio.com/reviews/submit/NBP3-38076C)

Earn gift cards/discounts by submitting a publication using this product:  
[www.novusbio.com/publications](http://www.novusbio.com/publications)

