

# Product Datasheet

## Proteasome 20S beta 3 Antibody [DyLight 550] NBP3-38629R

Unit Size: 0.1 ml

Store at 4C in the dark.

[www.novusbio.com](http://www.novusbio.com)



[technical@novusbio.com](mailto:technical@novusbio.com)

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:  
[www.novusbio.com/NBP3-38629R](http://www.novusbio.com/NBP3-38629R)

Updated 3/24/2025 v.20.1

Earn rewards for product  
reviews and publications.

Submit a publication at [www.novusbio.com/publications](http://www.novusbio.com/publications)

Submit a review at [www.novusbio.com/reviews/destination/NBP3-38629R](http://www.novusbio.com/reviews/destination/NBP3-38629R)



**NBP3-38629R**

Proteasome 20S beta 3 Antibody [DyLight 550]

Product Information	
Unit Size	0.1 ml
Concentration	Please see the vial label for concentration. If unlisted please contact technical services.
Storage	Store at 4C in the dark.
Clonality	Polyclonal
Preservative	0.05% Sodium Azide
Isotype	IgG
Conjugate	DyLight 550
Purity	Affinity purified
Buffer	50mM Sodium Borate
Product Description	
Host	Rabbit
Gene ID	5691
Gene Symbol	PSMB3
Species	Human, Mouse, Rat
Immunogen	<p>Recombinant fusion protein containing a sequence corresponding to amino acids 1-205 of human Proteasome 20S beta 3 (NP_002786.2).</p> <p>Sequence:  MSIMSYNGGAVMAMKGKNCVAIAADRRFGIQAQMVTTDFQKIFPMGDRLYIGL  AGLATDVQTV AQRLKFRLNLYELKEGRQIKPYTLMSMVANLLYEKRFGPYYTEP  VIAGLDPKTFKPFICSLDLIGCPMVTDDFVVS GTCAEQMYGMCESLWEPNMDP  DHLFETISQAMLNAVDRDAVSGMGVIVHIEKDKITTRTLKARMD</p>
Notes	<p>DyLight (R) is a trademark of Thermo Fisher Scientific Inc. and its subsidiaries.</p>
Product Application Details	
Applications	Western Blot, ELISA, Immunohistochemistry, Immunohistochemistry-Paraffin
Recommended Dilutions	Western Blot, ELISA, Immunohistochemistry, Immunohistochemistry-Paraffin
Application Notes	Optimal dilution of this antibody should be experimentally determined.





### **Novus Biologicals USA**

10730 E. Briarwood Avenue  
Centennial, CO 80112  
USA  
Phone: 303.730.1950  
Toll Free: 1.888.506.6887  
Fax: 303.730.1966  
nb-customerservice@bio-techne.com

### **Bio-Techne Canada**

21 Canmotor Ave  
Toronto, ON M8Z 4E6  
Canada  
Phone: 905.827.6400  
Toll Free: 855.668.8722  
Fax: 905.827.6402  
canada.inquires@bio-techne.com

### **Bio-Techne Ltd**

19 Barton Lane  
Abingdon Science Park  
Abingdon, OX14 3NB, United Kingdom  
Phone: (44) (0) 1235 529449  
Free Phone: 0800 37 34 15  
Fax: (44) (0) 1235 533420  
info.EMEA@bio-techne.com

### **General Contact Information**

www.novusbio.com  
Technical Support: nb-technical@bio-  
techne.com  
Orders: nb-customerservice@bio-techne.com  
General: novus@novusbio.com

### **Products Related to NBP3-38629R**

---

NBP2-24891R	Rabbit IgG Isotype Control [DyLight 550]
NBP1-72544-100ug	Recombinant Human Proteasome 20S beta 3 His Protein
NBL1-14878	Proteasome 20S beta 3 Overexpression Lysate
NB100-143	Nbs1 Antibody

---

### **Limitations**

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit [www.novusbio.com/guarantee](http://www.novusbio.com/guarantee)

Earn gift cards/discounts by submitting a review: [www.novusbio.com/reviews/submit/NBP3-38629R](http://www.novusbio.com/reviews/submit/NBP3-38629R)

Earn gift cards/discounts by submitting a publication using this product:  
[www.novusbio.com/publications](http://www.novusbio.com/publications)

